




- **Comfort** - Moisture & air permeable.
- **Working pressure** - Short stretch for a high working pressure.
- **Gap resistant** - Special notches and minimal overlying straps resist gapping and bulk.
- **EZ pull tab** - Aid for hassle free application.
- **Patent Pending 3 layer system** - Helps cool and wick moisture.

- Juzo Cooling Compression Wrap
- Juzo Half Circle Hook Fasteners
- Juzo Full Foot Liner

- Remove all the half circle hook fasteners from the straps.
- Machine wash warm.
- A garment bag is strongly recommended.
- Tumble dry low.
- **DO NOT USE:** bleach, chlorine, fabric softeners, or other laundry additives.
- Store garments in a dry place at room temperature.



1. The Juzo Cooling Compression Wrap comes with a Full Foot Wrap Liner.
 - If using the Juzo Wrap Liner, slide the liner over the foot and heel, pulling up the leg and over the calf. Position the top of the liner just below the knee, making sure the liner is smooth from top to bottom.
2. Attach the hook fasteners to the straps on the gray side of the Cooling Compression Wrap.
 - 2a. The straight end of the fastener should overlap the strap's end by 3/4 of inch. Press firmly to secure the straight end of the fastener to the strap.
 - 2b. The half-circle portion of the fastener should extend beyond the strap's end and attach to the body of the wrap.
3. Gently fold the fasteners over onto themselves to prevent the fasteners from inadvertently attaching to the garment. Each fastener can then be unfolded before attaching.
4. While in a resting position, i.e. sitting down, position the wrap on the leg just below the knee. (2 fingers below the knee). The seam should run down the center of the backside of the leg.
5. Secure the top strap of the wrap by counter pulling the top strap and the top EZ Pull at the same time. The strap should be just tight enough to gently hug the leg and hold the garment in place.
6. Continue working down the leg and fasten the remaining straps using enough tension for the straps to gently hug the leg. We recommend simultaneously counter pulling the two adjacent straps.
7. Once each strap is in place, start at the bottom of the leg and in unison pull the bottom strap and EZ Pull in opposite directions.
8. Continue working up the leg and counter pull the second and third straps. Make sure there are no gaps between straps. This can be best achieved by having adjacent straps slightly overlap and by inserting one strap into the notch.
9. Pull the top strap and the EZ Pull until it is comfortably tightened.
10. The wrap should have a snug fit and not be uncomfortable or binding. Throughout the day, if you feel any discomfort or the wrap feels too loose, then repeat steps 6 through 8, adjusting the level of tension that is applied.

Stay Comfortable with the Juzo Cooling Wrap