

Thank you for choosing Juzo, an international leader in medical compression garments since 1912.

Manufactured with care and detail, your Juzo Cooling Compression Wrap has several features that provide therapy, comfort, versatility and superior fit, including:


- **Comfort** - Moisture & air permeable.
- **Working pressure** - Short stretch for a high working pressure.
- **Gap resistant** - Special notches and minimal overlying straps resist gapping and bulk.
- **EZ pull tab** - Aid for hassle free application.
- **Patent Pending 3 layer system** - Helps cool and wick moisture

Included In This Package:

- Juzo Cooling Compression Wrap
- Juzo Hook Fasteners

Care and Maintenance

- Remove all the half circle hook fasteners from the straps.
- Machine wash warm.
- A garment bag is strongly recommended.
- Tumble dry low.
- **DO NOT USE:** bleach, chlorine, fabric softeners, or other laundry additives.
- Store garments in a dry place at room temperature.



Return Policy

Should it be necessary to return your garment for evaluation, Juzo requires that the garment be laundered (for hygienic reasons) prior to returning to the original place of purchase.

Fabric

49% Polyurethane / 45% Nylon / 6% Spandex • Not made with natural rubber latex.

Instructions for Use

It is recommended that you consult with your medical professional and fitter for proper guidelines on the use and application of this product. If you experience any discomfort, immediately remove the short-stretch compression wrap and consult with your medical professional..

Indications

Primary & secondary varices, chronic venous insufficiency, post-thrombotic syndrome, venous stasis ulcers, pre and post-vein procedure and lymphatic edema.

Contraindications

Absolute: peripheral arterial disease, decompensated heart diseases, septic phlebitis, phelgmasia coerulea dolens.

Relative: Leg infection, intolerance to compression or garment fabric, sensory disturbance to the limb.

Cautions

No side effects are currently known. However, if a negative result should occur (i.e. skin irritations) from the use of our products, please discontinue use of the garment and consult with a medical professional. The manufacturer is not liable for damage /injuries which have been caused by improper handling or misuse.

Steps to Putting on the Foot Wrap

1. Attach the hook fasteners to the straps on the gray side of the Cooling Compression Wrap.
 - 1a. The straight edge of the fastener should overlap the strap's end by 3/4 of inch. Press firmly to secure the straight end of the fastener to the strap.
 - 1b. The half-circle portion of the fastener should extend beyond the strap's end and attach to the body of the wrap.
2. Gently fold the fasteners over onto themselves to prevent the fasteners from inadvertently attaching to the garment. Each fastener can then be unfolded before attaching the strap to the garment.
3. Slide the wrap on your foot so that the heel fits into the heel cup. The seam at the back of the garment, above the heel opening should be centered at the back of the heel. If being worn with a Juzo Cooling Compression Wrap Calf Segment, it should be applied second and slightly overlap the bottom part of the calf segment.
4. Simultaneously pull the bottom strap, nearest the toes, and the bottom EZ Pull until the strap is comfortably snug.
5. After each strap is initially fastened, go back and pull the second and third straps again until they are both comfortably tightened around the heel and instep area.
6. Pull the top strap, nearest the ankle, and the top EZ Pull until it is comfortably tightened around the ankle.
7. The wrap should have a snug fit and not be uncomfortable or binding. Throughout the day, If you feel any discomfort or the wrap feels too loose, then repeat steps 4 through 6, adjusting the level of tension that is applied.

Note: If you are using a Juzo Cooling Compression Wrap Calf Segment, see care and application directions for use.

COOLING
COMPRESSION WRAP
Short Stretch • Comfortable

Venda de compresión de elasticidad reducida
Bandage compressif à extensibilité courte

- Moisture Wicking
- Cool To Touch
- Quick Dry

FOOT WRAP

brrr°
PROVEN TO KEEP YOU COOLER

INC00LFOOTWR-0525
les jambes en bonne santé et en mouvement.
Bandage compressif à extensibilité courte pour aider à améliorer la circulation et garder
Venda de compresión de elasticidad reducida para piernas saludables y siempre en movimiento.

Donning & Dörfling
Please Recycle
Printed in USA

For over 100 years, our family philosophy has been that medical garments should not inhibit a person's lifestyle.

Manufactured By
Julus Zorn, Inc
3690 Zorn Drive
Cuyahoga Falls, OH 44223

Distributed by Juzo Canada LTD.

See inside of insert for material
contents & proper care.

PCSZ-01-52A

Note: This size chart is not intended to replace a professional fitting by a qualified fitter.
The Juzo Compression Wrap will offer a varying degree of pressure that is dependent on the tension of the garment and the position and activity level of the wearer. Typical resting pressure will range between 20-40 mmHg and working pressure will range between 30-60 mmHg.

	(A1) Arch	(B) Ankle							
S	19-24 cm 7 1/2"- 9 1/2"	15-23 cm 6"- 9"	Reg	16 cm / 6 1/4"	18 cm / 7 1/4"	19.5 cm / 7 3/4"	Long	21 1/8 1/4"	19 cm / 7 1/2"
M	24-29 cm 9 1/2"-11 1/2"	20-28 cm 7 3/4"- 11"	Reg	18 cm / 7 1/4"	19.5 cm / 7 3/4"	Long	22.5 cm / 9"	23 cm / 9 1/4"	23 cm / 9 1/4"
L	29-34 cm 11 1/2"-13 1/4"	25-33 cm 9 3/4"-13"	Reg	19.5 cm / 7 3/4"	20.5 cm / 8 1/4"	Long	24 cm / 9 3/4"	25 cm / 9 3/4"	26 cm / 10"
XL	34-40 cm 13 1/4"-15 3/4"	30-40 cm 11 3/4"- 15 3/4"	Reg	20.5 cm / 8 1/4"	21.5 cm / 8 3/4"	Long	25.5 cm / 10"	26.5 cm / 10 1/4"	27.5 cm / 10 3/4"

Our patented 3-layer system features a Rapid Dry Outer Layer that wicks away moisture while the Ventilated Foam Middle Layer enhances airflow for a quick dry. The brrr° Fiber Inner Layer provides a cool feel against your skin.

Stay Comfortable with the Juzo Cooling Wrap

Moisture Wicking + Quick Dry + Cool To Touch

RAPID DRYING + VENTILATING FOAM + COOLING FIBER