ROHO®

Hybrid Select Cushion

Operation Manual

permobil

s h a p e f i t t i n g t e c h n o l o g y®

Contact Information

 ROHO, Inc.
1501 South 74th St.
Belleville, IL 62223 U.S.A.
U.S.A.:
800-736-0925
Fax 888-551-3449
Customer Support:
orders.roho@permobil.com
Outside the U.S.A.:
618-277-9150
Fax 618-277-6518
intlorders.roho@permobil.com
permobil.com

Warning

Indicates that not following the specified procedure could lead to potentially hazardous conditions resulting in serious injury.

Caution

Indicates that not following the specified procedure could lead to potentially hazardous conditions resulting in minor to moderate injury or damage to the equipment or other property.

Incident Reporting

In the event of a serious incident or injury, contact Customer Support.

Contents

Intended Use	2-3
Important Safety Information	
Parts Detail	4
ISO FLO Memory Control	
Offloading Option	
Setup	5-6
Cover and Cushion Components	7
Cleaning and Disinfecting	8
Troubleshooting	9
Product Specifications	
Storage, Transport, Disposal and Recycling	
Limited Warranty	
Español	10
Français	19
Deutsch	28
Italiano	37
Nederlands	46
Dansk	55
Svenska	64
Norsk	73
Suomi	82
Português	91
Hrvatski	100
Čeština	109
Polski	118
Slovenščina	127
Arabic	136
Hebrew	145

Symbols Glossary

	Manufacturer		FITS	Fits		SN	Serial number		FRONT	Front of product		CL	Cleaning
	Date of manufacture		REF	Model number			Weight limit		BACK	Back of product		DI	Disinfecting
	Consult instructions for use		QTY	Quantity		WDH	Size (width x depth x height)		BTM	Bottom of product		✓	Meets flammability
	LOT	Batch code		MD	Medical device		EC REP	Authorized representative in the European Community		CH REP	Authorized representative for Switzerland		

DO NOT use the ISO FLO Memory Control as a handle; the product may tear.

Please Note

Supplier: This manual must be given to the user of this product.

Operator (Individual or Caregiver): Before using this product, read the instructions and save for future reference.

Products may be covered by one or more U.S. and foreign patents and trademarks, including DRY FLOATATION®, ISO FLO® Memory Control, QUADTRO SELECT®, ROHO® and shape fitting technology®. Zytel® is a registered trademark of DuPont.

Product information is changed as needed; current product information is available at permobil.com. Contact Customer Support for previous versions of product information.

ROHO, Inc. is part of Permobil.

© 2020, 2022 Permobil

Rev. 2022-07-25

Intended Use Statement

The ROHO Hybrid Select Cushion (Hybrid Select) is an adjustable, air-filled, wheelchair support surface combining a ROHO DRY FLOATATION cellular design cushion insert with a contoured foam base, and is intended to conform to an individual's seated shape in the pelvic loading area to provide skin/soft tissue protection, positioning, and an environment to facilitate wound healing. Optionally, through removal of the IT Air Pad insert, the Hybrid Select is also intended for offloading of the ITs. There is a weight limit of 500 lbs. (227 kg), and the cushion must be properly sized to the individual. The Hybrid Select must be used with the supplied cover.

Contraindication: Based on clinical, scientific or engineering evidence, the Hybrid Select may not be suitable for individuals who require accommodation of pelvic asymmetry greater than 1 in. / 2.5 cm.

Compatibility: The Hybrid Select is - the ROHO Contour Base.
 NOT COMPATIBLE with the following: - the standing function available on some power wheelchairs.
 - 10 degrees or more of anterior tilt available on power wheelchairs.
 - more than 110 degrees of backrest recline on power wheelchairs or manual wheelchairs.

A clinician who is experienced in seating and positioning should determine whether the product is appropriate for the individual's particular seating needs. The product's intended use statement alone should not be relied upon to make that determination.

The medical products featured in this document are intended to be only one part of an overall care regimen that includes all seating and mobility equipment and therapeutic interventions. A clinician who is experienced in seating and positioning should identify this care regimen after evaluating the individual's physical needs and overall medical condition.

A clinician should also evaluate for visual, reading, and cognitive impairment to determine the need for caregiver assistance or the need for other assistive technology, such as large-print instructions, to ensure proper product use.

Important Safety Information

Warnings:

- Skin/soft tissue breakdown can occur due to a number of factors, which vary by individual. Check skin frequently, at least once a day. Redness, bruising, or darker areas (when compared to normal skin) may indicate superficial or deep tissue injury and should be addressed. If there is any discoloration to skin/soft tissue, STOP USE immediately. If the discoloration does not disappear within 30 minutes after disuse, immediately consult a healthcare professional.
- Check inflation frequently, at least once a day.
- DO NOT use a product that is underinflated or overinflated, because 1) the product benefits will be reduced or eliminated, resulting in an increased risk to skin and other soft tissue, and 2) the individual may become unstable and vulnerable to falling. Carefully follow all "Setup" instructions. If the product does not appear to be holding air, or if you are not able to inflate or deflate the product, see "Troubleshooting". Immediately contact the equipment provider, distributor, or Customer Support if the problem persists.
- DO NOT use the components of this product separately. Assemble and use this product only as instructed in this manual.
- The cushions, foam and cover MUST be compatible sizes and MUST be used as directed in this manual. IF THEY ARE NOT:
 - 1) The benefits provided by the product may be reduced or eliminated, increasing the risk to skin or to other soft tissue, and
 - 2) the individual may become unstable and vulnerable to falling.

Warnings:

- DO NOT use the cushion on top of, or in conjunction with, any other products or materials, except as indicated in this manual. Using non-compatible products or materials, whether inside or outside the cover: 1) may reduce or eliminate the benefits of the cushion, increasing the risk to skin or to other soft tissue, and 2) may cause the individual to become unstable and vulnerable to falling.
- DO NOT place obstructions between the individual and the cushion. The cushion and the cover MUST be compatible sizes and MUST be used as directed in this manual. Except for the compatible cover provided with this product, placement of any items between the individual and the cushion: 1) may reduce or eliminate the benefits of the cushion, increasing the risk to skin or to other soft tissue, and 2) may cause the individual to become unstable and vulnerable to falling.
- DO NOT use any part of the product as a water flotation device (e.g. a Life Preserver). It will NOT support you in water.
- The cushion surface will acclimate to the surrounding temperature. Take appropriate precautions, especially when the cushion will contact unprotected skin.
- DO NOT expose the product to high heat, open flames, or hot ashes. Testing or certification claims, including for flammability, may no longer apply to this device when it is combined with other products or materials. Check testing and certification claims for all products used in combination with this device.
- If you are unable to perform any tasks described in this manual, seek assistance by contacting the clinician, equipment provider, distributor, or Customer Support.

Important Safety Information, continued**⚠ Warnings:**

- DO NOT use the Hybrid Select with a power wheelchair at greater than 10 degrees of anterior tilt. The product benefits will be reduced or eliminated, resulting in an increased risk to skin and other soft tissue.
- DO NOT recline greater than 110 degrees in any wheelchair while seated on a Hybrid Select. The product benefits will be reduced or eliminated, resulting in an increased risk to skin and other soft tissue.
- When reclining in any wheelchair while seated on a Hybrid Select: after reaching the desired recline angle, check to make sure there is no pinching or tissue compression.
- Areas of contact with the Hybrid Select may shift when the wheelchair is adjusted to a recline angle, which may result in an increased risk of tissue injury. BOTH after reaching the desired recline angle AND after returning to an upright position: 1) perform a Transition Check and Immersion Check on the Hybrid Select, as instructed in this manual; 2) if using the Hybrid Select without the IT Air Pad (Offload Option), check the position of the ITs to confirm desired offloading.
- A Hybrid Select that is not securely attached to a power wheelchair may shift during wheelchair adjustments, resulting in a safety risk to the user. Make sure that the Hybrid Select is securely installed on the wheelchair. Confirm that the number and position of hook and loop fasteners provide secure attachment to the seat pan.

⚠ Cautions:

- Do not use a pump, cover or repair kit other than those compatible with the product.
- DO NOT modify the cushion or any components. Doing so may result in damage to the product and will void the warranty.
- DO NOT allow the air cells to come into contact with oil-based lotions or lanolin, which may degrade the material.
- Prolonged exposure to ozone may degrade materials used in the cushion, affect the performance of the cushion, and void the product warranty.
- Changes in altitude may require adjustment to the cushion.
- Keep the cushion away from sharp objects.
- If the cushion has been in temperatures lower than 32°F (0°C) and exhibits unusual stiffness, allow the cushion to warm to room temperature.
- Periodically, check the cover and components for wear and damage and replace if necessary.
- DO NOT use any of the following as a handle for carrying or pulling the cushion: inflation valve, hoses, or the ISOFLO Memory Control. Carry the cushion by the cushion base or by the carrying handle on the cover. Using any other components to carry or pull the cushion may damage the product.

⚠ Warnings - Motor Vehicle Transportation:



- Failure to pay attention to these warnings could result in severe injury to the individual in the wheelchair or to others.
- Whenever possible, transfer out of the wheelchair and into a manufacturer-installed vehicle seat without the cushion, and use the vehicle's crash-tested occupant restraint system.
- The cushion must be correctly and securely installed, used, and maintained according to the instructions and warnings provided by the cushion manufacturer.
- When the wheelchair is used as the vehicle seat, the cushion should only be used with wheelchair frames that are equipped with securement points specified in ANSI/RESNA WC-4:2017, section 19 (WC19) or ISO 7176-19, and that have been successfully crash tested on a wheelchair in the forward-facing position.
- The wheelchair with the cushion must be used with an effective wheelchair securement system and a properly positioned, crash-tested pelvic and shoulder-belt restraint, or Wheelchair Tiedown and Occupant Restraint System (WTORS), following the manufacturer's instructions.
- If the cushion has been involved in an accident during transport, check the cushion for damage and replace if necessary. See the "Troubleshooting" section in this manual.
- For more information about safely transporting wheelchair users, refer to: <http://wctransportation-safety.umtri.umich.edu/ridesafe-brochure>.
- A ROHO QUADTRO SELECT Cushion was included in a wheelchair seating system that was dynamically tested by a third-party testing facility for use in a motor vehicle. The system tested meets all applicable criteria for wheelchair seating systems proposed in Section 5.1 of ANSI/RESNA WC-4:2017 Section 20 "Wheelchair Seating Systems for Use in Motor Vehicles."
- Test results apply to seating systems intended for use by adults and children with a total body mass greater than 23 kg (51 lb, or approximately 6 years of age) who are transported facing forward in their wheelchairs in all types of motor vehicles.
- To determine whether motor vehicle third-party testing has been performed with a specific seating system combination, refer to all provided equipment operation manuals or contact the equipment manufacturer.

Parts Detail


Package contents: Hybrid Select Cushion, Hybrid Select Cover, hand pump, repair kit, operation manual, limited warranty supplement, product registration card.

A serial number label is located on the bottom of the air pads and the foam base.

Note: To order replacement parts, contact the product provider or Customer Support.

ISO FLO Memory Control


Offloading Option

⚠ IMPORTANT! If you are using the Hybrid Select without the IT Air Pad (as an offloading cushion) follow these instructions, and then continue with the Cushion Setup instructions. Consult with a clinician to adjust the cushion as prescribed.

1. Remove the cover. See the "Cover and Cushion Components" section in this manual.

2. Locate the tubing that attaches the IT Air Pad to the ISO FLO, and detach the quick disconnect.

3. Pull the tubing through the opening in the foam base and remove the IT Air Pad.

4. Make sure the Large Air Pad is positioned correctly on the foam base, with the hook/loop strips carefully aligned.

5. Store the IT Air Pad in a container that will protect it from moisture, contaminants, wear and damage. See additional storage instructions in this manual.

Setup

STOP HERE if using the Hybrid Select without the IT Air Pad. Refer to the "Offloading Option" in this manual before completing the Setup instructions on this page.

⚠ Warnings - Before Use: Read all safety information and perform all steps below. Transfer to the cushion only after performing all steps included in "Setup" (two pages). Make sure the cushion is properly sized to the individual.

SETUP - PART ONE: Remove the cover and unlock the ISO FLO.

- Using a stable surface, such as a table, unzip and open the cover. Turn the product over and make sure the quick disconnects are securely connected to the ISO FLO.
Note: When attaching the quick disconnects: The tubes can be attached in any order.


Inflation valve

Quick Disconnects
ISO FLO Memory Control

- Make sure the ISO FLO is unlocked. See the "ISO FLO Memory Control" section.

Note: When the ISO FLO is unlocked, the air flows freely through all of the air cells. You will not be able to inflate the cushion unless the ISO FLO is unlocked.

SETUP - PART TWO: Perform atmospheric inflation.

⚠ IMPORTANT! Make sure there is NO WEIGHT on the air pads. If the individual is sitting on the cushion, transfer off.

- Add air:
 - Open the inflation valve.
Notes: To open, turn at least one full rotation. You can use the hand pump nozzle to open the inflation valve.
 - Inflate with the hand pump until all of the air cells feel firm.
Note: All air cells will inflate at the same time.
- Adjust to atmospheric pressure:
 - Remove the hand pump, and DO NOT close the inflation valve.
 - Allow air to escape from the inflation valve. Listen (or place a hand in front of the inflation valve) to see if the airflow has stopped.
- Close the inflation valve completely. The cushion is now inflated to atmospheric pressure.

SETUP - PART THREE: Cover and place in the wheelchair.

- Before closing the cover:
 - Make sure the Large Air Pad is positioned correctly, with the hook and loop strips aligned.
 - Make sure the Large Air Pad is securely attached to the foam base.
 - Make sure there is NO GAP between the Large Air Pad and the foam base.
- Zip the cover closed (see cover installation instructions).

No gap!

- Place the Hybrid Select on the wheelchair (air cells facing up and non-skid surface of the cover facing down).

Note: Use the symbol on the cover to help orient the cushion.

⚠ Warning: To ensure the product performs as intended: Make sure that the back corners of the cushion contact the back canes of the wheelchair. For wheelchairs that do not have back canes, make sure the cushion is flush with the backrest upholstery. Check for interference with other wheelchair components. Make sure that the wheelchair foot supports and back support are correctly positioned.


Continue with the instructions on the next page.

Setup, continued

SETUP - PART FOUR: Transferring to the wheelchair (and POSITIONING, if needed):

Before transferring, make sure the inflation valve is closed. The individual should transfer to the wheelchair and sit in their usual sitting position.

Note: The inflation valve should be in the front, on the left side of the seated individual.

If a clinician determines that POSITIONING is needed: After inflating all air cells to atmospheric pressure and transferring to the wheelchair: Make sure the ISO FLO is unlocked. Perform the prescribed positioning. Then, lock the ISO FLO to maintain the position.

See the ISO FLO Memory Control section in this manual.

SETUP - PART FIVE: Transition and Immersion Checks


⚠ Warning: Skin/soft tissue injury may result if the transition from air cells to the foam base is not level, due to overinflated or underinflated air cells. The air cells must not be overinflated. Refer to the warning about overinflation/underinflation in the "Important Safety Information" section of this manual.

⚠ IMPORTANT! It may be necessary to re-check and re-adjust the Transition and Immersion areas, to ensure proper inflation.

Reminder! The ISO FLO must be open while adding or removing air.

FIRST: Transition Check: Make sure the transition from the air cells to the foam is level, in the area indicated.


Instructions: Place a hand beneath the individual, where the air cells meet the foam base (indicated by the oval shape), and feel the transition from air pad to foam. Follow the instructions below.

Proper Inflation

It feels LEVEL where the air cells meet the foam. There is no bump or ridge of overlapping air cells and foam.


What to do: Perform an immersion check, as shown below.

Too much air in the air cells

It feels like a bump or ridge of air cells lying on top of foam. The air cells are too high.

What to do: Remove a small amount of air and check again. Repeat as needed until the air cells are properly inflated. Then, perform an immersion check, as shown below.


Not enough air in the air cells

It feels like a bump or ridge of foam. The air cells are not level with the foam.


What to do: Add a small amount of air and check again. Repeat as needed until the air cells are properly inflated. Then, perform an immersion check, as shown below.

SECOND: Immersion Check: After the air cells are properly inflated to ensure a level Transition from air cells to foam, check to make sure there is enough air (or space) under the ITs.

With the IT Air Pad: Place a hand beneath the lowest bony prominences (ITs) and make sure there is at least ½" of air. Add or remove a small amount of air, if needed. Repeat until the air cells are properly inflated.

Proper Inflation
(with the IT Air Pad)

Underinflated
(with the IT Air Pad)

Without the IT Air Pad (Offloading Option): Place a hand beneath the lowest bony prominences (ITs) and make sure that the ITs do not contact the foam. Add or remove a small amount of air, if needed. Repeat until the air cells are properly inflated.

Proper Inflation
(without the IT Air Pad)

Underinflated
(without the IT Air Pad)

Cover and Cushion Components

To remove the cover:


Unzip the cover. Make sure the inflation valve is inside the cover. Push the ISOFLO into the slot and inside the cover. Place the product, air cells facing up, on a stable surface. Remove the cover.

To remove the cushion components:

1. Detach all quick disconnects from the ISOFLO.

2. Pull the tubing through the opening in the foam base.

3. Remove the cushion components, releasing them from the hook/loop fastener strips.

To install the cushion components:

1. If using the IT Air Pad, place it into the cutout area of the foam base. The tubing should face the opening in the foam base.


2. Place the Large Air Pad onto the foam base, carefully matching the hook/loop strips. Make sure the Large Air Pad is securely attached to the foam base. Make sure there is NO GAP between the Large Air Pad and the foam base.

⚠ IMPORTANT! The flange of the IT Air Pad (if using) should lie under the flange of the Large Air Pad.

The "flange" is indicated by the dotted line.

3. Insert the tubing through the opening in the foam base.
4. Using the quick disconnects, attach each tube to the ISOFLO.

When attaching the quick disconnects: The tubes can be attached in any order.

5. Position the inflation valve tubing in the channel.

To install the cover:


1. Confirm that the cushion components are installed correctly and connected to the ISOFLO.

2. Place the foam base on a stable surface. Slide the cover over the front of the cushion, and make sure it is snug.

3. Pull the ISOFLO out through the front slot in the cover. Push the inflation valve through the hole in the corner of the cover. Make sure the tubing is in the channel.

4. Place the cover around the rear corners of the foam base. Close the zipper.

⚠ IMPORTANT! Check to make sure the cover fits correctly on the contours of the cushion.

Note: If used, hook fastener strips that were removed prior to washing should be replaced.

Cleaning and Disinfecting

⚠ Warnings:

- Cleaning and disinfecting are separate processes. Cleaning must precede disinfection. If the product becomes soiled, or before use by a different individual: clean, disinfect, and check the product for proper functioning.
- Clean the product regularly and when it becomes soiled.
- Follow all manufacturer bleach container and germicidal disinfectant safety guidelines.
- Note about sterilization: High temperatures accelerate aging and will damage the product assembly. The products featured in this manual are not packaged in a sterile condition, nor are they intended or required to be sterilized prior to use. If institution protocol demands sterilization: First, follow the cleaning and disinfecting instructions; then, open the inflation valve and use the lowest possible sterilizing temperature, but not to exceed 70°C (158°F), for the shortest time possible. DO NOT autoclave.

⚠ IMPORTANT! After cleaning and disinfecting and before use, follow all instructions to set up the cushion.

Recommendation: When washing by hand, use room-temperature water.

Cover

Before cleaning and disinfecting, remove the cover from the cushion. Remove any hook fastener pieces from the bottom of the cover. Close the cover flap.

To clean the cover: Machine wash cold with mild detergent, gentle cycle; or wipe clean with neutral detergent and water. Tumble dry low or air dry.

To disinfect the cover: Hand wash, using 1 part household liquid bleach per 9 parts water. Keep the product wet with the bleach solution for 10 minutes, and then rinse thoroughly. Or machine wash in hot water (60°C). Tumble dry low or air dry.

Cushion and Foam Components

⚠ Caution: DO NOT machine wash or machine dry the cushions or foam.

Before cleaning and disinfecting: Remove the cover. Remove the cushion components from the foam base. (See the "Cover and Cushion Components" section in this manual.)

Large Air Pad and IT Pad

⚠ Cautions: Use of the following may cause damage to the cushion components: abrasives (steel wool, scouring pads); caustic, automatic dishwasher detergents; cleaning products that contain petroleum or organic solvents, including acetone, toluene, Methyl Ethyl Ketone (MEK), naphtha, dry cleaning fluid, adhesive remover; oil-based lotions, lanolin; ozone gas; ultraviolet light cleaning methods. DO NOT expose the cushion or foam components to direct sunlight. Rinse thoroughly; detergent residue may cause the air cells to stick together.

⚠ Warning: DO NOT allow water or cleaning solution to enter the cushion components. DO NOT submerge the tubing in any liquid, including water.

To clean the cushion components: Hand wash, using mild liquid hand soap, hand dishwashing liquid, laundry detergent, or a multipurpose detergent (follow product label instructions.) Gently scrub all cushion surfaces, using a soft plastic bristle brush, a sponge or a cloth. Rinse with clean water. Air dry.

To disinfect the cushion components: Repeat cleaning instructions, using 1 part household liquid bleach per 9 parts water. Keep the product wet with the bleach solution for 10 minutes. Rinse thoroughly with clean water. Air dry.

⚠ Caution: Rinse thoroughly. Detergent residue may cause air cells to stick together.

Foam Base

⚠ Warning: DO NOT submerge the foam in any liquid, including water.

⚠ Warning: Disinfectants are not effective on porous surfaces such as foam. If the foam becomes soiled or contaminated, DO NOT use with multiple individuals. Replace soiled or contaminated foam.

To clean the foam: To clean, wipe the foam gently with a cloth dampened with water and a household detergent, then wipe again with a clean, slightly damp cloth. Allow the foam to air dry thoroughly before installing the cover.

Care Symbols

Hand wash.

Wipe with a damp cloth.

Machine wash gentle, cold.

Do not machine dry.

Tumble dry low.

Line dry in the shade.

Machine wash hot, normal, at the temperature shown.

Do not machine wash.

Bleach (1 part bleach: 9 parts water).

Dry flat in the shade.

Troubleshooting

For additional assistance, contact the equipment provider, distributor, or Customer Support.

Not holding air	» Confirm that the inflation valve is completely closed. » Inspect all components for damage. Look for holes in the cushion. If very small holes or no holes are visible, follow the instructions in the Repair Kit provided with the product. For damage to the ISOFLO or the tubing or for large holes or leaks in the cushion components, see the Limited Warranty supplement, or contact Customer Support. »  Repair Kit Warning: DO NOT allow any liquid to enter the cushion components.
Not able to inflate or deflate the cushion	» Confirm that the ISOFLO is open. » Confirm that the inflation valve is open. » Make sure that all quick disconnects are securely connected.
Uncomfortable / unstable	» Confirm that the air pad has the correct amount of air. See the Transition Check and Immersion Check instructions. » Make sure the cushion components are correctly placed on the foam base; the hook / loop fastener strips on the air cells and the foam base should match up. » Make sure the air cells face up and the top of the cover faces up. » Confirm that the cushion is centered beneath the individual and that the cushion is properly sized to the individual and the wheelchair. If the cushion still feels uncomfortable or unstable after checking all of the above, immediately consult the prescribing healthcare professional.
The cushion slides on the wheelchair	» Make sure the non-skid bottom of the cushion cover faces down. » Use directional labels on the cushion and on the cover to orient the cushion. » When possible, use hook and loop fasteners to further secure the cushion to the seating surface. Hook fasteners are available upon request.
A component is damaged	Contact Customer Support for a replacement.

Product Specifications

Materials: Large Air Pad / IT Pad: neoprene rubber; made without latex; Inflation valve: Zytel; Tubing: urethane; ISOFLO Memory Control: polycarbonate, aluminum, polypropylene; Foam Base: polyurethane; Cover: two-way stretch, fluid-resistant top and sides, with a zipper on the back and sides, and a non-skid bottom. Refer to cover care label for material content and for compliance to flammability requirements.

WIDTH (W) - in. (cm) ➡	14 (35.5)	15 (38)	16 (40.5)	17 (43)	18 (45.5)	20 (51)	Notes: - All measurements are approximate. - Width and depth measurements are taken from the top of the foam base and are rounded to the nearest ¼ in. or 0.5 cm. - Item numbers and cover labeling indicate the wheelchair size that the product is designed to fit. - ** Depths shown include ¾ in. (1.5 cm) of foam base that fits between the wheelchair back posts. WEIGHT: 4 lbs. / 1.8 kg (based on 1818) SEATED HEIGHT: 3½ in. / 9 cm
DEPTH (D)** - in. (cm) ⬇	Contoured Foam Base width x depth						
14¾ (37.5)	1414						
15¾ (40)		1515	1615				
16¾ (42.5)	1416	1516	1616	1716	1816	2016	
17¾ (45)		1517	1617	1717	1817	2017	
18¾ (47.5)			1618	1718	1818	2018	
20¾ (52.5)			1620	1720	1820	2020	

Cushion expected lifetime: 5 years

Cover expected lifetime: 1 year

Accessories* Compatible with Hybrid Select: ROHO Planar Solid Seat Insert

* Not available in all markets. Contact Customer Support.

Storage, Transport, Disposal and Recycling

Storage or Transport: Clean and disinfect the product. Open the inflation valve. Place in a container that protects the product from moisture, contaminants and damage. After storing or transporting, inspect the product to ensure it was not damaged (see "Troubleshooting"), and follow instructions to prepare the product for use.

Disposal: The components of the products in this manual are associated with no known environmental hazards when properly used and when disposed of in accordance with all local/regional regulations. At the end of its useful life, treat as healthcare waste and dispose of according to local healthcare-waste guidelines. Incineration must be performed by a qualified, licensed waste management facility.

Recycling: Contact a local recycling agency to determine recycling options for the product.

Limited Warranty

Limited Warranty Term from the date the product was originally purchased: cushion, 36 months; cover, 6 months. The warranty does not apply to punctures, tears or burns. See also the Limited Warranty supplement provided with the product, or contact Customer Support.

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Permobil LTD