



Products / Wound Care / Bandages / Conforming Stretch Gauze / Conforming Stretch Gauze, 2" x 4.1 yds



Conforming Stretch Gauze, 2" x 4.1 yds 602

12/BG, 8 BG/CS

Product Features

Rayon/poly knitted stretch gauze bandage

Provides slight compression to the site while securely holding to any body contour

Specifically designed not to slip while in use

Unique weave construction and finished edges prevent unraveling and linting

Attributes

Brand Name	Dukal
Size	2" x 4.1 yd
HCPC Code	A6442
Sterile	No
Weight	1.45 lbs

Packaging

Packaging	12/BG, 8 BG/CS
Inner Packaging	12/BG
Case Packaging	8 BG/CS
Total Case Quantity	96
Inner Packaging Barcode	(01)10665973015108
Case Packaging Barcode	(01)20665973015105



About Dukal

[Who We Are](#)

[Our Legacy](#)

Solutions

[Innovation](#)

[Supply Chain](#)

[QA/RA](#)

[Contract Manufacturing](#)