



Premium Care New Sponges, 4-Ply, 2" x 2"

6112

200/BG, 40 BG/CS

Product Features

- Rayon/poly blend construction
- The effective combination of weight and open aperture allows each sponge to wick quicker and offer superior absorbency
- Virtually non-linting

Attributes

Brand Name	Dukal
Size	2" x 2"
Ply	4-Ply
HCPC Code	A6216
Sterile	No
Weight	8.2 lbs

Packaging

Packaging	200/BG, 40 BG/CS
Inner Packaging	200/BG
Case Packaging	40 BG/CS
Total Case Quantity	8000
Inner Packaging Barcode	(01)10665973014026
Case Packaging Barcode	(01)20665973014023

