

Compreflex Standard

ARM

Compression Level

20–30mmHg



Colors



Black

Beige

Features

The Compreflex Standard Arm is an inelastic compression device for the arm. The device features an attached inner sleeve for easy application and the patented 90° Flexi-Zone™ for improved movement at the elbow. The Compreflex Standard Arm may be paired with the Dorsal Pocket Glove or Secure Gauntlet for compression in the hand.

Ideal for


- Lymphedema
- Post-surgical edema
- Non-specific edema

Includes

- Compreflex Standard Arm
- 1 Pair of non-compressive Cotton Liners
- 1 Butterfly Strap

Product Selection Tip

May be right for patients with lymphedema of the arm



Note: Hand compression wear sold separately



Dorsal Pocket Glove



Secure Gauntlet

Compreflex Standard

ARM

Patient Name:	Contact Name:
Account Name:	Account #:
Account Phone #:	P.O. #:
Ship Name & Address:	Date:

Product Information

Product includes one Compreflex Standard Arm and one pair of Cotton Liners.

<input type="checkbox"/> Right Arm	Size:	Item #:
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length:	Quantity:
<input type="checkbox"/> Left Arm	Size:	Item #:
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length:	Quantity:

Companion Products (sold separately)

<input type="checkbox"/> Dorsal Pocket Glove	Size:	Item #:
<input type="checkbox"/> Right Hand	Quantity:	
<input type="checkbox"/> Dorsal Pocket Glove	Size:	Item #:
<input type="checkbox"/> Left Hand	Quantity:	

Dorsal Pocket Glove Size Chart (A- E)

(Available in Black Only)

	Small	Medium	Large	X Large
E	15-19cm	19-22cm	22-26cm	26-29cm
A	14-20cm	16-22cm	18-24cm	20-26cm
Black				
Right	1201-DGR-R	1202-DGR-R	1203-DGR-R	1204-DGR-R
Left	1201-DGR-L	1202-DGR-L	1203-DGR-L	1204-DGR-L

Circumference

Left Right


C	_____	_____
B1	_____	_____
B	_____	_____
A	_____	_____
E	_____	_____

D Length

Short: 38-42.9cm
Regular: 43-47.9cm
Long: 48-53cm

Measure length along
dorsal aspect

D	_____	_____
---	-------	-------



Compreflex Standard Arm Size Chart (A-C)

	Small	Medium	Large	X Large
C	23-32cm	28-37cm	33-43cm	39-49cm
B1	22-30cm	26-34cm	30-39cm	35-44cm
B	20-27cm	24-31cm	28-35cm	32-39cm
A	14-19cm	15-20cm	16-21cm	18-23cm
Black				
Short	1401-ARS	1402-ARS	1403-ARS	1404-ARS
Regular	1401-ARR	1402-ARR	1403-ARR	1404-ARR
Long	1401-ARL	1402-ARL	1403-ARL	1404-ARL
Beige				
Short	1411-ARS	1412-ARS	1413-ARS	1414-ARS
Regular	1411-ARR	1412-ARR	1413-ARR	1414-ARR
Long	1411-ARL	1412-ARL	1413-ARL	1414-ARL



Coolflex Standard ARM

Compression Level

20–40mmHg



Color



Black

Features

The Coolflex Standard Arm is ideal for individuals managing lymphedema following mastectomy or other surgical intervention. The garment is made of a lightweight, cooling material that transfers heat and releases moisture away from the skin, leaving it cool and dry. The Coolflex material has two layers which gives it a cushioned feel and helps to protect and provide support the limb.

The product features an easy to grip closure system with hooks for simple, independent donning. The unique pleated elbow design contains a stretch panel to encourage a full range of motion for comfortable all-day wear.

Pair the Coolflex Standard Arm with the Dorsal Pocket Glove for compression to the hand.

Ideal for

- Lymphedema
- Post-surgical edema
- Non-specific edema

Includes

- Coolflex Standard Arm
- 1 Pair of non-compressive Cotton Liners

Product Selection Tip

Select the Coolflex Standard Arm for post mastectomy patients



Note: Hand compression wear sold separately



Dorsal Pocket Glove



Secure Gauntlet

Coolflex Standard

ARM

Patient Name:	Contact Name:
Account Name:	Account #:
Account Phone #:	P.O. #:
Ship Name & Address:	Date:

Product Information

Product includes one Coolflex Standard Arm and one pair of Cotton Liners.

<input type="checkbox"/> Right Arm	Size:	Item #:
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length:	Quantity:
<input type="checkbox"/> Left Arm	Size:	Item #:
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length:	Quantity:

Companion Products (sold separately)

<input type="checkbox"/> Dorsal Pocket Glove	Size:	Item #:
<input type="checkbox"/> Right Hand	Quantity:	
<input type="checkbox"/> Dorsal Pocket Glove	Size:	Item #:
<input type="checkbox"/> Left Hand	Quantity:	

Dorsal Pocket Glove Size Chart (A- E)

(Available in Black Only)

	Small	Medium	Large	X Large
E	15-19cm	19-22cm	22-26cm	26-29cm
A	14-20cm	16-22cm	18-24cm	20-26cm
Black				
Right	1201-DGR-R	1202-DGR-R	1203-DGR-R	1204-DGR-R
Left	1201-DGR-L	1202-DGR-L	1203-DGR-L	1204-DGR-L

Circumference

Left Right


C	_____	_____
B1	_____	_____
B	_____	_____
A	_____	_____
E	_____	_____

D Length

Short: 38-42.9cm
Regular: 43-47.9cm
Long: 48-53cm

Measure length along dorsal aspect

D	_____	_____
---	-------	-------



Coolflex Standard Arm Size Chart (A-C)

(Available in Black Only)

	Small	Medium	Large	X Large
C	23-32cm	28-37cm	33-43cm	39-49cm
B1	22-30cm	26-34cm	30-39cm	35-44cm
B	20-27cm	24-31cm	28-35cm	32-39cm
A	14-19cm	15-20cm	16-21cm	18-23cm

Right

Short	2101-ARS-R	2102-ARS-R	2103-ARS-R	2104-ARS-R
Regular	2101-ARR-R	2102-ARR-R	2103-ARR-R	2104-ARR-R
Long	2101-ARL-R	2102-ARL-R	2103-ARL-R	2104-ARL-R

Left

Short	2101-ARS-L	2102-ARS-L	2103-ARS-L	2104-ARS-L
Regular	2101-ARR-L	2102-ARR-L	2103-ARR-L	2104-ARR-L
Long	2101-ARL-L	2102-ARL-L	2103-ARL-L	2104-ARL-L

