

Airgo Comfort-Plus Lightweight Transport Chair

HCPCS:



- New and improved with more security and comfort enhancing features
- Flip-back detachable arms make it easier to transfer in and out, or simply get closer to a table
- Comfort-Grip 8" wheels provide shock absorption and better brake performance
- Anatomically correct seat-to-armrest height for a natural posture
- Fold-down back for easy transportation and storage
- Removable swing-away footrests with tool-free adjustment of length for greater comfort
- Heel-loops and seatbelt included as standard features for increased safety
- Dual locking rear wheels for safe, secure transfers
- Chip-resistant finishing and no sharp edges on hardware
- Lightweight and compact

PRODUCTS	Description
700-841	Airgo 17" TC, Black, 1 EA,2CV,Ret.
700-846	Airgo 19" TC, Black, 1 EA,2CV,Ret.
700-848	AIRGO 19" ALUM TC PLAID, 1EA, 2CV, Ret

SPECIFICATIONS	Description
Warranty	Frame and Cross-brace 3 Years, All Other Metal Parts 6 Months, Upholstery, Plastic and Rubber 3 Months
Product Dimensions	26" (L) x 22" 24" (W) x 37" (H)
Seat Dimensions	17" 19" (W) x 16" (D)
Handle Height	37"
Weight Capacity	300 lb
Seat To Floor Height	19"
Folded Dimensions	26" (L) x 10.5" (W) x 37" (H)
Weight	20.5 lb