

Bed

Hälsa & Hälsa Plus

Features at a glance

- Standard and Plus sizes available
- Trendelenburg & Reverse Trendelenburg
- Dual posts on each side increases stability at all heights
- Low position: 9.875"
- High position: 25.625"
- Disassembles into four section for easy transport or storage
- Standard size: 36" x 80"; Plus size: four width options - 36", 39", 42", and 48"
- Weight capacity of up to 375 lb for standard and 500 lb for Plus/bariatric model



Bed

Hälsa & Hälsa Plus

Adjustable electric home comfort bed

The Hälsa Bed is an electric, adjustable home care bed available in standard and Plus sizes and includes a 10-button hand pendant. The Trendelenburg (and Reverse Trendelenburg) feature comes standard and allows clients to lay flat on their back on a 15–30 degree incline with their feet elevated above the head. Additionally, it offers the capability of moving both the head and foot sections simultaneously, an upgraded large castor option, and easy disassembly for transportation and storage.



Standard features:

- Natural or wine colour head/footboard
- Hooks for incontinence bag
- 3 locking castors & 1 directional locking castor
- Side and foot mattress retainers
- Receptacle for trapeze bar
- 2 receptacles for IV holder stems
- Virtually noise-free with low electrical consumption DC Linak motors

Optional features:

- Assist bar or half side rails, interchangeable between head and foot sections
- Full side rails (only fits standard bed)
- New! Slip-on rail pad covers (shown)
- Trapeze bar
- Battery backup
- Upgraded large castor option



Rail pad