


Parts List

Joerns[®] Lifting & Repositioning **Hoyer[®] Advance[®]**

To avoid injury, read user's manual before using.



| Item | Part # | Description | Qty |
|------|--------------|--|-----|
| 1 | 0Y0066 | 100mm rear caster (braked) | 2 |
| 2 | 0Y0109 | 75mm front caster (non-braked) | 2 |
| 3 | 0Y0033 | Mast extrusion (inc. labels) | 1 |
| 4 | 0Y0034 | Boom extrusion (inc. boom end caps and labels, except product name decal) | 1 |
| 5 | 0Y0035 | Push handle assembly (inc. fixings) | 1 |
| 6 | 0Y0036 | Casted boom to mast pivot joint (inc. connection bolts/fixings/spacers) | 1 |
| 7 | 0Y0037 | Spreader bar assembly complete (6-point) | 1 |
| 8 | 0Y0309 | Spreader bar / cradle quick release pin kit | 1 |
| 9 | 0Y0379 | Spreader bar washer and o-ring spares kit | 1 |
| 10 | 0Y0038 | Base assembly (inc. base casting, bottom plates, welded post & labels, not internal parts - excluding fixings) | 1 |
| 11 | 0Y0039 | Left leg extrusion (inc. fixings) | 1 |
| 12 | 0Y0040 | Right leg extrusion (inc. fixings) | 1 |
| 13 | 0Y0106 | Control box - Standard | 1 |
| 14 | 400-3003 | Battery pack | 1 |
| 15 | 0Y0107 | User handset (Standard control box) | 1 |
| 16 | 0Y0031 | Actuator - 350mm stroke (inc. location pins, sleeves, washers, circlips and labels) | 1 |
| 17 | 0Y0069 | Battery and controller bracket (inc. fixings) | 1 |
| 18 | 0Y0041 | Foot pedal mechanism complete | 1 |
| 19 | 0Y0042 | Foot push pad assembly complete | 2 |
| 20 | 0Y0108 | Boom to mast pivot plastic caps | 2 |
| 21 | 0Y0045 | Front caster leg inserts including fixings (plastic) | 2 |
| 22 | 0Y0046 | Knob kit (inc. mast tightening knob, elastic leg restraint and leg disengage handle) | 1 |
| 23 | 0Y0047 | Fixings kit (all bolts, screws and fixings) | 1 |
| 24 | 0Y0149 | Actuator fixing kit | 1 |
| 25 | 0Y0729 | Boom end plastic cap | 1 |
| 26 | 0Y0032 | Hydraulic pump - 350mm stroke | |
| 27 | 400-3016 | Cable Power (CNTRL BOX) Linkak HME | 1 |
| - | 0Y0044 | Decal set (exc. serial number labels) (not shown) | 1 |
| - | 0Y0070 | Handset clip (not shown) | 2 |
| - | 808003.10013 | Sling retainer assembly (not shown) | 2 |

