

Bed

Tröst & Tröst Plus



Features at a glance

- Standard and Plus sizes available
- Trendelenburg, Reverse Trendelenburg, and Comfort chair positioning
- 12-button hand pendant
- Low position: 6.75"; High position: 30"
- Flip down directional castor locks at head end, plus swivel locking castor
- Centre step safety lock at foot end
- Mattress retainers and wall bumpers included
- Trost STD and Plus T/RT is $+17^{\circ} / -17^{\circ} \pm 3^{\circ}$
- Standard size: 36" x 80"; Plus size: 39" x 80" or 42" x 80" and length extension to 84"
- Weight capacity of up to 600 lb



Hand control

Bed

Tröst & Tröst Plus

Adjustable long-term care bed

The Tröst Bed is an electric, adjustable bed ideal for the long-term care setting and community. It offers a low height of 6.75" and a high height of 30" to meet all your caregiving and client needs. Tröst Bed has a 12-button hand pendant and optional footboard attendant controls, so customized comfort is an easy adjustment. The Trendelenburg (and Reverse Trendelenburg) feature comes standard and allows clients to lay flat on their back on a 15–30 degree incline with their feet elevated above the head. Available in a natural or wine colour head/footboard (*sold separately, except when bought as a package*), that are easy to use with quick release hardware.



Setup video

Optional features:

- Plus model expandable to 39" and 42"
- Length extension from 80" to 84" (added at foot end only)
- Embedded footboard attendant controls
- Half side rails, interchangeable between head and foot sections
- Tall half side rails
- Assist bar, interchangeable between head and foot sections
- Trapeze bar
- **New!** Slip-on rail pad cover
- Battery backup
- Transport trolley



Wheel-locks



New! Slip-on rail cover