

Tröst

LITE Bed



Features at a glance

- Standard size of 36" x 80" with optional 4" length extension (80" - 84")
- 8-button hand pendant
- Low position: 7.3"; High position: 28.3"
- Four locking casters, two of which are flip down directional caster locks at the head section
- Mattress retainer and wall bumper included
- Maximum Back Angle: 65°±3°, Maximum Knee Angle: 30°±3°
- Weight capacity: Up to 500 lb
- Virtually noise-free, smooth operation
- Two-piece bed frame for easy delivery and tight spaces



Tröst

LITE Bed

Adjustable long-term care bed

The Tröst LITE is a 3 function, electric, adjustable bed, ideal for the long-term care setting and community. With a two-piece low-profile bedframe, a lower bed position of 7.3" and a max height of 28.3" this long-term care bed will fill all your caregiving needs. An 8-button hand pendant and optional staff controls allow for easy and customizable comfort. Available with a natural or wine colour head/footboard (sold separately, except when bought as a package.)



Optional features:

- Length extension from 80" to 84" (added at foot section only)
- Headboard and footboard (optional colors)
- Staff Control (embedded, hanging, or both)
- Standard or tall half side rails (head section, foot section, or both)
- Assist bar (head section, foot section, or both)
- Full side rails
- New! Slip-on rail pad covers (shown)
- Battery backup
- Transport trolley

The **Tröst Lite Bed (HSR)** package includes: 36" head/footboard and a set of half side rails, while the **Tröst Lite Bed (FSR)** package includes: 36" head/footboard and a set of full side rails.

Hand control

