

Tubigrip®

Multi-purpose tubular support bandage

Provides firm support for sprains, strains and swelling



What size to choose?

Product code	Tubigrip size	Minimum limb circumference (cm)	Minimum limb circumference (inch)	Maximum limb circumference (cm)	Maximum limb circumference (inch)	Sub-bandage pressure (mmHg)
1435	A	10	3.9	14.5	5.7	20 - 30
1436	B	14.5	5.7	29	11.4	10 - 20
1443	C	22	8.7	35	13.8	
1437	D	26	10.2	40	15.7	
1434	E	29	11.4	46	18.1	
1438	F	35	13.8	53	20.9	
1439	G	40	15.7	62	24.4	5 - 15
1440	J	54	21.3	82	32.3	

Tubigrip is available under a full range of sizes to suit most applications. The sub-bandage pressures are obtained in double layer.

Reference: Mölnlycke Health Care. Data on file.

Find out more at www.molnlycke.ca

Mölnlycke Health Care Canada. 2010 Winston Park Drive, Suite 100 Oakville, ON L6H 5R7. Phone: 905-829-1502
© 2023 Mölnlycke Health Care AB. All rights reserved. HQIM002303.

