

AeroWalk™ Ultra-Lite Rollator

Elevate Every Step with the AeroWalk™ Ultra-Lite Rollator



Refined Style, Designed to Move with You
Elegant colors & woodgrain accents highlight the seamless, weld-free frame.



Ultra-Lightweight
Less than 15 pounds*
with a compact
side-to-side fold



Optimal Support
Extendable, adjustable-
height backrest provides
lumbar support at rest



Smooth Mobility
10" front wheels
& 8" back wheels for
indoor-outdoor use



Customized Comfort
Height-adjustable
handles & ergonomic
hand grips for all-day use



Convenient Storage
Stylish, removable
bag keeps
essentials close

Lightweight mobility meets modern design with the AeroWalk™ Ultra-Lite Rollator. Combining a seamless, weld-free frame and innovative back support, it sets the new benchmark for portable performance at **less than 15 pounds***

Item #	Description	UPC	UOM
10277BL	AeroWalk Ultra-Lite Rollator, Tall, Blue	822383027104	1/ea
10277GY	AeroWalk Ultra-Lite Rollator, Tall, Gray	822383027098	1/ea
10277WT	AeroWalk Ultra-Lite Rollator, Tall, White	822383027111	1/ea
10277BL-S	AeroWalk Ultra-Lite Rollator, Standard, Blue	822383031590	1/ea
10277GY-S	AeroWalk Ultra-Lite Rollator, Standard, Gray	822383031606	1/ea
10277WT-S	AeroWalk Ultra-Lite Rollator, Standard, White	822383031613	1/ea

Product Specs

Material	Aluminum	Product Weight	14.8 lb.*
Wheel Size	10" front 8" back	Shipping Weight	18.5 lb.
Handle Height	31.5" - 35.4" Standard 31.5" - 37.25" Tall	Weight Capacity	330 lb.
Unit Dimensions	23.25" (W) x 26.5" (D)	HCPC	E0143 + E0156
Seat Dimensions	17.75" (L) x 10" (W) x 23" or 20" (H)	Warranty	Limited Lifetime
Seat-to-Floor Height	20" Standard 23" Tall		

*Without storage bag.



Easy Side-to-Side Folding
Allows for
effortless travel
and storage