

## AeroCarbon Ultra-Lite Rollator

HCPCS: E0143 + E0156, E0143 + E0156



- Effortless Maneuverability, Travel & Storage: Easy to lift, carry, and maneuver at just 11.5 pounds. Stores compactly in closets, trunks, and more thanks to a slim fold.
- Durable & Easy to Maintain: Provides structural integrity without the bulk with carbon fiber construction. Reduces common fracture points with a weld-free frame and supports simplified maintenance with fewer replacement parts.
- Style Customers Love: Converts even hesitant mobility aid users with a sleek design and sophisticated matte finish.
- Confidence & Comfort: Offers stable performance indoors and outdoors with 9.5" front wheels and 8" rear wheels, plus snag-free use with enclosed brake cables. Customize support with adjustable handle height and 17.5" wide seat with ergonomic backrest.

### PRODUCTS

### Description

10278BK	Aerocarbon Rollator, Black, 1/ea
10278GY	Aerocarbon Rollator, Grey, 1/ea

### SPECIFICATIONS

### Description

Product Weight	11.5 lb
Material	Carbon Fiber
Seat Dimensions	17.5" x 9"
Handle Height	31.7"37.6"
Weight Capacity	330 lb
Height	31.7"
Length	26.77"
Seat To Floor Height	23"
Casters	9.5"
Overall Width	23.62"
Width Between Arms	17.72"