

Airgo Transport Chair

HCPCS:



- Lightweight and compact, almost 50% lighter than conventional wheelchairs
- Fold-down back for storage
- Black nylon upholstery, with padded armrests
- 8" mag-style wheels with dual locking rear wheels
- Swing-away removable footrests, and standard safety belt
- Easy to fold and convenient to store

PRODUCTS	Description
700-850	Airgo Steel TC, 17" Wide, Chrome, 1 EA
700-855	Airgo Steel TC, 19" Wide, Chrome, 1 EA
SPECIFICATIONS	Description
Warranty	3 Year Limited
Product Dimensions	23" (L) x 21", 23" (W) x 39.7" (H)
Seat Dimensions	17" or 19" (W) x 16" (D)
Handle Height	39.7"
Weight Capacity	300 lb (136 Kg)
Carton Shipping Weight	28.6 lb
Casters	8"
Folded Dimensions	23" (L) x 11.4" (W) x 39.7" (H)
Weight	28.6 lb net