

AeroWalk Ultra-Lite Rollator

HCPCS: E0143 + E0156,E0143 + E0156



- Ultra-Lightweight Design with sleek aluminum frame weighs less than 15 lbs.
- Seamless frame has no visible welds and three sleek color options to choose from.
- Concealed brake cable for improved durability and safety.
- Adjustable handle height 31.5" - 35.4" (Standard) and 31.5" - 37.3" (Tall).
- Backrest is adjustable for premium lumbar support
- Large 10" front wheels provide optimal stability and rolling comfort.
- Durable and comfortable seat for rest on the go.
- Easy side-to-side folding for storage and transport.
- Removable storage bag to easily carry personal items.

PRODUCTS	Description
10277BL	Aerowalk Rollator, Blue, 1/ea
10277BL-S	Aerowalk Rollator, Standard, Blue, 1/ea
10277GY	Aerowalk Rollator, Gray, 1/ea
10277GY-S	Aerowalk Rollator, Standard, Gray, 1/ea
10277WT	Aerowalk Rollator, White, 1/ea
10277WT-S	Aerowalk Rollator, Standard, White, 1/ea

SPECIFICATIONS	Description
Product Weight	14.8 lb (without storage bag)
Material	Aluminum
Seat Dimensions	17.7" (W) x 10" (L)
Handle Height	Standard: 31.5-35.4 Tall: 31.5" - 37.3"
Weight Capacity	330 lb
Height	31.5"
Length	26.6"
Seat To Floor Height	Standard: 20" Tall: 23"
Casters	10"
Overall Width	23.3"
Width Between Arms	19.3"