

Customer and Order Information

Purchase Order *	Tag	v.0101
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		* required

Base Chair

Fuze T50n 50° Manual Tilt Wheelchair

- | | | | |
|-------------------------------------|---------|--|----------|
| <input checked="" type="checkbox"/> | 2601000 | Fuze T50 50° Manual Tilt Wheelchair | \$ 4,300 |
| <input checked="" type="checkbox"/> | 2601008 | Narrow Frame Upgrade Reduces overall chair width to seat width plus only 4". | \$ 200 |

Frame Type & Weight Capacity

Fixed Width Frame Select a Fixed Seat Width below.

- | | | | |
|--------------------------|---------|--|--------|
| <input type="checkbox"/> | 2602010 | 300 lbs Weight Capacity | |
| <input type="checkbox"/> | 2602011 | 400 lbs Weight Capacity Not compatible with Transit Tie Down System. | |
| <input type="checkbox"/> | 2602025 | Adjustable Width Frame 250 lbs Weight Capacity. Select an Adjustable Seat Width below. | \$ 100 |

Frame Options

- | | | | |
|--------------------------|---------|---|--------|
| <input type="checkbox"/> | 2602222 | Take-Apart Frame The Take-Apart frame allows attendants to take the Fuze apart in seconds without tools. pdgmobility.com/take-apart | \$ 500 |
|--------------------------|---------|---|--------|
- This feature reduces the weight of the heaviest piece of the wheelchair to 29 lbs or less. Perfect for transporting it in the trunk of a car or SUV. Not compatible with Power Tilt, Ventilator Tray, Ventilator Battery Tray, DRB, User Wheel Locks. Take-Apart rear crossbar is not width adjustable. PDG will provide one no-charge crossbar in a different width over the lifetime of each chair ordered with adjustable width frame (can be ordered with the chair or later). Note: Crossbar for the T50n is 4" shorter than the Seat Width.



- | | | | |
|--------------------------|---------|---------------------------------------|--|
| <input type="checkbox"/> | 2602300 | Transit Tie Down System (TTDS) | |
|--------------------------|---------|---------------------------------------|--|
- TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.

Seat Width

Overall chair width is +4.0" for 12" RW

Fixed Seat Widths Order with Fixed Frame.

- | | | | | | | | | |
|--------------------------|---------|----------------|--------------------------|---------|----------------|--------------------------|---------|----------------|
| <input type="checkbox"/> | 2603016 | 16" Seat Width | <input type="checkbox"/> | 2603017 | 17" Seat Width | <input type="checkbox"/> | 2603018 | 18" Seat Width |
| <input type="checkbox"/> | 2603019 | 19" Seat Width | <input type="checkbox"/> | 2603020 | 20" Seat Width | <input type="checkbox"/> | 2603021 | 21" Seat Width |
| <input type="checkbox"/> | 2603022 | 22" Seat Width | | | | | | |
| <input type="checkbox"/> | 2603023 | 23" Seat Width | <input type="checkbox"/> | 2603024 | 24" Seat Width | | | \$ 375 |

Adjustable Seat Widths Order with Adjustable Frame.

- | | | | | | |
|--------------------------|---------|----------------------------|--------------------------|---------|----------------------------|
| <input type="checkbox"/> | 2603216 | 16"-19" Range - Set at 16" | <input type="checkbox"/> | 2603217 | 16"-19" Range - Set at 17" |
| <input type="checkbox"/> | 2603218 | 16"-19" Range - Set at 18" | <input type="checkbox"/> | 2603219 | 16"-19" Range - Set at 19" |
| <input type="checkbox"/> | 2603319 | 19"-20" Range - Set at 19" | <input type="checkbox"/> | 2603320 | 19"-20" Range - Set at 20" |

Seat Depth

Standard Seat Depths (Adjustable)

☐ 2604015 15" Seat Depth

☐ 2604016 16" Seat Depth

☐ 2604017 17" Seat Depth

☐ 2604018 18" Seat Depth

☐ 2604019 19" Seat Depth

☐ 2604020 20" Seat Depth

Custom Seat Depths

Requires Stroller Handle Back - Straight. Not compatible with Reclining or Fold-down Back.

\$ 400

☐ 2604021 21" Seat Depth

☐ 2604022 22" Seat Depth

Frame Colour

PDG Premium Colours

View colour swatches on PDG's website. pdgmobility.com/colours

☐ 2605001 Black Shimmer

☐ 2605015 Sky Blue

☐ 2605016 Poppy Red

☐ 2605003 Fine Texture Black

☐ 2605010 Lollypop Blue

☐ 2605011 Lollypop Red

☐ 2605017 Arctic Grey

☐ 2605018 Tropical Teal

☐ 2605019 Tangerine Orange

☐ 2605020 Olive Green

☐ 2605021 Lush Lilac

☐ 2605012 Bubble Gum Pink

☐ 2605022 Leaf Green

☐ 2605023 Boysenberry Purple

☐ 2605014 Daffodil Yellow

☐ 2605200 Custom Colour

Contact PDG for specifications.

\$ 300

Armrests

☐ 2606000 No Armrests

Height Adjustable Armrests (Single Post)


\$ 295

☐ 2606114 10"-14" Range

☐ 2606116 12"-16" Range


☐ 2606100 Fixed Height Armrests (12")

Full Length Armpad only.

 Armrest Height

If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.

Seat Pan to top of Armrest: _____ Does not include height of Armpad

 Armrest Position

Indicate preference for Armrest position. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated

☐ 2606402 Mid

☐ 2606403 Front

Locking Flip-up Armrests

Not recommended for users with excessive strength or tone. Not compatible with Fold-down Back.

\$ 295

☐ 2606310 Set at 10"

☐ 2606311 Set at 11"

☐ 2606312 Set at 12"

☐ 2606313 Set at 13"

☐ 2606314 Set at 14"

Armpads

Regular Armpads

☐ 2607100 Desk Length - 2" x 10"

☐ 2607101 Full Length - 2" x 14"

Waterfall Armpads

\$ 115

☐ 2607105 Desk Length - 2" x 10"

☐ 2607106 Full Length - 2" x 14"

Wide Armpads

May interfere with wheels in some tilt configurations.

\$ 200

☐ 2607110 Desk Length - 3" x 10"

☐ 2607111 Full Length - 3" x 12"

Gel Armpads

May interfere with wheels in some tilt configurations.

\$ 135

☐ 2607120 Full Length - 2" x 12"

☐ 2607121 Full Length - 3.5" x 12"

Contour Arm Trough with Hand Pad

\$ 125

Not compatible with Lap Tray, Locking Flip-up Armrests. May interfere with wheels in some tilt configurations.

☐ 2607130 Right (4.5" x 18")

☐ 2607131 Left (4.5" x 18")

Armrest Options

Armrest Wideners

Move Armrests outward by up to 1" on each side. Seat size does not change.

☐ 2607210 1" Total Additional Width

\$ 95

☐ 2607211 2" Total Additional Width

\$ 195

Seat Pan

☐ 2608000

No Seat Pan

☐ 2608001

Solid Aluminum Seat Pan

\$ 135

Seat Height

➤ Don't forget

Select desired Seat-to-Floor Height

STF Height is adjustable.

12" Rear Wheels

Caster size	16"	17"	18"	19"	20"
5"	<input type="checkbox"/>	<input type="checkbox"/>			
6"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
8"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
6" Soft Roll	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Casters

Urethane Casters

☐ 2611051

5" x 1.25" Caster size

☐ 2611061

6" x 1.25" Caster size

☐ 2611081

8" x 1.25" Caster size

\$ 85

☐ 2611063

6" x 1.50" Caster size

\$ 135

☐ 2611084

8" x 2.00" Caster size

\$ 135

Pneumatic Casters

☐ 2611161

6" x 1.25" Caster size

\$ 110

☐ 2611181

8" x 1.25" Caster size

\$ 120

☐ 2611184

8" x 2.00" Caster size

\$ 150

☐ 2611260


Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters

\$ 275

Rear Wheels

☒ 2612012

12" Composite Wheels



Rear Wheel Position

Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated

☐ 2612401

Rear

More stable

☐ 2612402

Mid

☐ 2612403

Front

Less stable

Axle Type

☐ 2613000

Quick Release Rear Axles

☐ 2613010

Fixed Rear Axles

Rear Tires

☐ 2614000

Urethane Tires

☐ 2614010

Pneumatic Tires

May result in a finished Seat Height 1/2" greater than stated.

\$ 60

☐ 2614020

Pneumatic Tires with Flat-Free Inserts

May result in a finished Seat Height 1/2" greater than stated.

\$ 95

Handrims

☒ 2615000

No Handrims

Wheel Locks & Brakes

User Wheel Locks

Not compatible with Take-Apart Frame, Foot Operated Wheel Lock, Rear-mounted Wheel Position.

\$ 150

☐ 2616010

Push to Lock

☐ 2616020

Pull to Lock

☐ 2616030

Attendant Wheel Locks (Push to Lock)

☐ 2616005

Foot Operated Wheel Lock

Rear-mounted unilateral single-activation lock. Not compatible with Ventilator Tray, User Wheel Locks.

\$ 595

☐ 2616300

Drum Brakes

Available with Urethane Tires only. Compatible with 17.5", 18.5", 19.5" STF Height only.

\$ 795

User-activated Unilateral Wheel Lock

Pull to Lock. Available up to 22" Seat Width only.

\$ 795

☐ 2616310

Left Activation

☐ 2616311

Right Activation

Wheel Lock Options

☐

2616100 Wheel Lock Extensions

\$ 60

Anti-Tippers

☒

2617010 Rear Anti-Tippers

Required Safety Option.

\$ 75

Back Posts

Back Post Style

☐ 2618100 Stroller Handle Back

☐ 2618105 Cane Handle Back

Not compatible with Headrest Mounting Plate, Power Tilt Manual Override.

Back Post Type

☐ 2618200 Straight

☐ 2618205 8° Bend

Back Post Height

Measured from the Seat Pan to the Handle.

☐ 2618420 Short (17")

☐ 2618425 Tall (23")

Back Post Folding

Fold-down Back is not compatible with Dynamic Rocker Back Interface, Reclining Backrest.

☐ 2619220 Fold-down

☐ 2619221 Non-folding

Seat to Back Angle

☐ 2619090 Set at 90°

☐ 2619096 Set at 98°

☐ 2619102 Set at 106°

☐ 2619108 Set at 114°

☐ 2619114 Set at 122°

Back Options

Angle/Height Adjustable Stroller Handle

Not compatible with Headrest Mounting Plate or Cane Handle Back.

\$ 150

☐ 2619010 Button on Outside

☐ 2619011 Button on Inside

☐ 2619050 Removable Cane Handle Push Bar

Not compatible with Stroller Handle Back.

\$ 140

Endomorph Back

Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12".

Not compatible with Flip-up Armrests.

\$ 395

☐ 2618102 1" Offset per side

☐ 2618103 1.5" Offset per side

☐ 2618104 2" Offset per side

☐ 2619212 Headrest Mounting Plate

Not compatible with Angle/Height Adjustable Stroller Handle, Cane Handle Back. Not required if mounting to a solid Back.

\$ 95

☐ 2619240 Seating Dynamics Dynamic Rocker Back Interface

\$ 2,995

Not compatible with Take-Apart Frame, Fold-down Back, 30° Recline, Standard Upholstery. Only availabe with Stroller Handle - 8° Bend.

Finished Back Post Height is 7" taller than selected Back Post Height. Maximim compative Seat Depth is 20". 90° and 98° Seat Angle only.

Back Upholstery

☐ 2620000 No Back Upholstery

Standard Upholstery

Not available with Cane Handle Back.

\$ 195

☐ 2620418 Short (18" Height)

☐ 2620421 Regular (21" Height)

☐ 2620450 Tension Adjustable Upholstery

Height of Upholstery is 16". Not available with Cane Handle Back.

\$ 295

Position of Upholstery

Select the desired position of the top of the Upholstery relative to the Seat Pan.

17" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>								

23" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Regular					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Tension Adj				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Tilt Activation

Attendant Tilt only

Select trigger mounting location below.

☐ 2623011 Right Back Post

☐ 2623012 Left Back Post

See next page for more options >

Power Tilt 40° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt

- | | | | | |
|--------------------------|---------|--|--|----------|
| <input type="checkbox"/> | 2623400 | Standard Power Tilt (250 lbs Weight Capacity) | | \$ 1,795 |
| <input type="checkbox"/> | 2623401 | Standard Power Tilt (350 lbs Weight Capacity) | Requires Fixed Frame. | \$ 2,095 |
| <input type="checkbox"/> | 2623410 | Power Tilt with Manual Tilt Override (250 lbs Weight Capacity) | Allows for instant manual release of tilt in case of power tilt failure. | \$ 2,295 |
| <input type="checkbox"/> | 2623411 | Power Tilt with Manual Tilt Override (350 lbs Weight Capacity) | Allows for instant manual release of tilt in case of power tilt failure. Requires Fixed Frame. | \$ 2,595 |
- Power Tilt Manual Override not compatible with Cane Handle Back.

- | | | | | |
|--------------------------------|---------------------------|---|---|--------|
| ➤ Select Control Switch | Toggle Switch | <input type="checkbox"/> 2623421 Installed on Right side | <input type="checkbox"/> 2623422 Installed on Left side | |
| | Button Switch | <input type="checkbox"/> 2623431 Installed on Right side | <input type="checkbox"/> 2623432 Installed on Left side | |
| | <i>Stealth</i> Egg Switch | <input type="checkbox"/> 2623440 Not installed | | |
| | Optional for Egg Switch: | <input type="checkbox"/> 2623445 Swing Away Mounting Hardware | Not compatible with all Headrest makes and models. | \$ 625 |
- Specify Headrest to be used: _____

Tilt Limit Optional. If no tilt limit is selected, the full tilt range will be available.

- | | | |
|---|---|---|
| <input type="checkbox"/> 2623220 Set at 20° | <input type="checkbox"/> 2623230 Set at 30° | <input type="checkbox"/> 2623240 Set at 40° |
|---|---|---|

Recline

- | | | |
|----------------------------------|------------------------------------|--------|
| <input type="checkbox"/> 2623310 | Reclining Backrest (Manual) | \$ 895 |
|----------------------------------|------------------------------------|--------|
- 30° maximum recline. Not compatible with Cane Handle Back, Fold-down Back. Not available with 22" Seat Depth.
Finished Back Post Height is 3" taller than selected Back Post Height.

Front Rigging

- | | | | | |
|----------------------------------|-------------------------|----------------------------------|---------------------------|--|
| <input type="checkbox"/> 2625000 | No Front Rigging | <input type="checkbox"/> 2625001 | Sector Blocks only | |
|----------------------------------|-------------------------|----------------------------------|---------------------------|--|
- Easy-Out™ Front Rigging** Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out
- | | | | |
|----------------------------------|--------------------------------------|---|--------|
| <input type="checkbox"/> 2625410 | Fixed 70° | | \$ 95 |
| <input type="checkbox"/> 2625420 | Universal Angle Adjustable (0°-110°) | Set at angle: _____ (use increments of 10°) | \$ 295 |
- Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.
- Swing Away Front Rigging** Legacy Design. Not compatible with Easy-Out. Requires pre-tilt of 5°.
- | | | | | |
|----------------------------------|----------------------------------|---|----------------|--------|
| <input type="checkbox"/> 2625130 | Swing Away 60° | <input type="checkbox"/> 2625140 | Swing Away 70° | |
| <input type="checkbox"/> 2625160 | Swing Away 90° | Includes 11"-16" Extension Tubes. Not compatible with 8" Casters. | | \$ 75 |
| <input type="checkbox"/> 2625200 | Swing Away Elevating Legrests | Includes Calf Pads. | | \$ 375 |
| <input type="checkbox"/> 2625300 | 1" Offset Hangers for Swing Away | | | \$ 75 |
- ☐ 2625110 **Standard Plug In 70° Front Rigging** Legacy Design. Not compatible with Easy-Out.

Footplates

- | | | | |
|----------------------------------|---|--|--------|
| <input type="checkbox"/> 2626010 | Composite Footplates | | |
| <input type="checkbox"/> 2626020 | Aluminum Angle Adjustable Footplates | \$ 165 | |
| <input type="checkbox"/> 2626030 | One Piece Fixed Footplate | Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable. | \$ 350 |
| <input type="checkbox"/> 2626040 | One Piece Flip-up Footplate | Requires selection of Composite or Aluminum Angle Adjustable. | \$ 395 |
| <input type="checkbox"/> 2626050 | Interlocking Angle Adjustable Footplates | Available for Seat Widths up to 21" only. | \$ 195 |

Footrest/Legrest Accessories

- | | | | |
|--|--|---|--------|
| Calf Pad Not compatible with Gel Padding. | | | |
| <input type="checkbox"/> 2627012 | For Easy-Out™ Front Rigging | Includes reinforced premium Calf Pad mount. | \$ 200 |
| <input type="checkbox"/> 2627010 | For Swing Away Front Rigging | | \$ 135 |
| <input type="checkbox"/> 2627020 | Gel Padding | Not compatible with Calf Pad or Easy-Out Front Rigging. | \$ 125 |
| Gel Knee Pad | | | |
| <input type="checkbox"/> 2627032 | For Easy-Out™ Angle Adjustable Front Rigging | | \$ 200 |
| <input type="checkbox"/> 2627030 | For Swing Away ELR Front Rigging | | \$ 100 |
| <input type="checkbox"/> 2627040 | Length Adjustable Heel Loops | | \$ 50 |
| <input type="checkbox"/> 2627050 | Bodypoint Padded Calf Strap | Sizes greater than 22" are not padded. | \$ 195 |

Front Rigging Extension Tubes

Regular Extension Tubes

No selection required for 90° Front Rigging

☐ 2628010

Short

14"-18" for Easy-Out Front Rigging

13"-17" for Swing Away and Standard Plug In Front Rigging


☐ 2628020

Long

18"-24" for Easy-Out Front Rigging

15"-19" for Swing Away and Standard Plug In Front Rigging

\$ 75



Lower Leg Length

If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.

Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.

Miscellaneous Options

Stealth Flip Down Abductor

Requires Solid Seat Pan.

\$ 650

☐ 2630060

Small Pad

3"x4"x3"x2" size

☐ 2630062

Large Pad

5"x5"x4"x3" size

☐ 2630070

QR Bracket - Accessory Mount

\$ 200

☐ 2630080

Gel Hip Guides

Reduces available Seat Width by 2".

\$ 495

☐ 2630090

Hip Belt Mount

Provides improved mounting options for Bodypoint Hip Belt or other third-party belts.

\$ 55

☐ 2630150

Ventilator Tray

Reduces ground clearance. Chair may be unstable when unoccupied. Minimum Seat Width is 20".
Not compatible with Power Tilt, Foot Operated Wheel Locks.

\$ 295

☐ 2630152

Ventilator Bracket

Compatible with LTV1 (laptop) ventilator.

\$ 295

☐ 2630160

Clear Lap Tray (Slide On)

Available for Seat Widths of 17"-22" only. Not compatible with Contour Arm Trough.

\$ 490

☐ 2630200

Q'Straint Lap Belt

Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length.

\$ 245

☐ 2630222

Wheelchair Backpack

Black with PDG Logo. 17"x14" size.

\$ 145

☐ 2630170

Adhesive Back Hook Strip, 2" x 30"

For Cushion mounting.

\$ 10

On-Chair Seating & Positioning



Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.