

Manual Wheelchairs

TiLite CR1


permobil






Key features

- Exceptionally lightweight frame
- ± 1" Center of Gravity fine tuning adjustment
- Integrated, fold-down backrest for easy transport
- Individually made, just for you



Performance components

- Specially designed caster bearings for optimal maneuverability and weight savings
- PBO fiber and titanium reinforcements
- New carbon fiber caster wheels



Custom profiled sideguard

- Supports backrest without need for backrest hinges
- Individually designed for tailored fit



Dual-angled front end

- Designed to increase environmental access
- Optimized caster placement



[12.0096, 12.0161]





Component colors

Anodized options



COG bracket & trim rings



Specifications

Transport weight	9 lb
Seat width	13" – 18"
Seat depth	13" – 19"
Front seat height	17.5" – 21"
Seat slope	7°, 10°, 13°
Front frame angle	85°
Seat back height	10" – 17"

Camber	0°, 2°, 4°
Center of gravity	2" – 5"
Backrest angle	80° - 101°
Compatible with SmartDrive	Yes
Weight limit	231 lb
Warranty	Limited