

TiLite X Canadian Order Form

Folding Frame

Effective June 2025 (New items in red)

Page 1 of 8

PROVIDER INFORMATION		RTS INFORMATION		
Dealer:	Account #:	RTS:	Phone:	
PO:	Ordered By:	Email:		
Email:		SHIP TO INFORMATION		
Mark For:	User Weight:	Name:	Ship Via:	
TiLite Sales Rep:		Address:		
TiLite Information:	Promo Code:	City:	State:	Zip:
Ship Date:	Order #:	Phone:	Fax:	
Please contact the Permobil Manual Customer Service Team at chairs.tilite@permobil.com - Phone: 800.736.0925 or Fax: 866.586.2413				

FRAME STYLE

- ☐ Standard Aluminum Frame (Full Adjustability, 300 Pound (136 Kg) Weight Limit) \$2795
- ☐ Standard Aluminum Frame with Transit Tie-Down (Full Adjustability, 300 Pound (136 Kg) Weight Limit) \$2995
- ☐ Heavy Duty Aluminum Frame (350 Pound (158 Kg) Weight Limit. Available with 18" to 22" Seat Widths and 16" to 20" Seat Depths)
Requires Aluminum Casters – Provided as NCO) \$995

SMARTDRIVE POWER ASSIST SYSTEM



- ☐ Permobil Smart Drive MX2+* (N/A with rear seat heights less than 13.5") \$7,989
 - ☐ SwitchControl Buttons NCO
 - ☐ SwitchControl Button with Mono Jack Option NCO
 - ☐ Buddy Button NCO
 - ☐ Carrying Bag \$110
 - ☐ Universal Padded Handle \$31
- *SmartDrive apps available on Apple App Store and Google Play Store for Apple Watch Series 5 or newer and Samsung Galaxy Watch 4 or newer.

TRANSIT OPTIONS

- ☐ Transit Tie-Down (Consists of four Tie-Down Brackets. N/A with HD Frame, Fixed Height Non-Removable Armrest or Swing Away Tubular Armrest.
(Please refer to the user manual for important information regarding the Transit Tie-Down Option) \$375
- ☐ Transport Bracket Option (Consists of four Tie-Down Brackets. Unoccupied Use Only. N/A with HD Frame, Fixed Height Non-Removable
Armrest or Swing Away Tubular Armrest) \$395

FRAME FINISHES

- | | | |
|---|---|--|
| <input type="checkbox"/> High Gloss Black
<input type="checkbox"/> Ocean Blue Metallic
<input type="checkbox"/> Midnight Blue Pearl
<input type="checkbox"/> Titanium Grey
<input style="color: red;" type="checkbox"/> Indigo Purple
<input type="checkbox"/> OD Green
<input style="color: red;" type="checkbox"/> Black Galaxy | <input style="color: red;" type="checkbox"/> Red Pepper
<input type="checkbox"/> Tangerine Metallic
<input style="color: red;" type="checkbox"/> White Pearl
<input style="color: red;" type="checkbox"/> Pomegranate Red
<input style="color: red;" type="checkbox"/> Hyper Blue
<input style="color: red;" type="checkbox"/> Silver Streak
<input style="color: red;" type="checkbox"/> Seafoam Green | <input type="checkbox"/> Neon Pink
<input type="checkbox"/> Acid Green Metallic
<input type="checkbox"/> Black Lava
<input style="color: red;" type="checkbox"/> Blue Mist
<input type="checkbox"/> Flame Orange
<input style="color: red;" type="checkbox"/> Lavender Dusk
<input style="color: red;" type="checkbox"/> Daffodil Yellow STD |
|---|---|--|

FRAME MEASUREMENTS (Please see information on page 8 regarding frame measurements)

For odd size seat depths, frame depth will be 1" longer than seat sling. Seat upholstery is 1" less than seat depth and begins 1" after front edge of back post.

Seat Width A										
		14	15	16	17	18	19	20	21	22
Seat Depth B	14	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
	15	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
	16	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	17	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	18	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	19		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	20			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Note: Configurations including Seat Widths 19"-22" \$365
 Configurations including 19" Seat Depth \$275



Permobil.com
 2701 West Court Street | Pasco, Washington 99301
 509.352.2245 Tel | 866.586.2413 Toll-Free
chairs.tilite@permobil.com

FRONT SEAT HEIGHT RANGES

Front Wheel Size	Front Seat Heights c
3" Front Wheel	13.5" -18"
4" Front Wheel	14.5" - 18.5"
5" Front Wheel	15.5" - 19"
6" Front Wheel	16.5" - 19.5"
8" Front Wheel	18.5" - 20.5"

FRONT SEAT HEIGHT c Measure from floor to top of seat tube at front edge of seat sling.

☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20" ☐ 20.5"**STD**

REAR SEAT HEIGHT RANGES

Rear Wheel Size	Rear Seat Heights d
22" (501)	13" - 18.5"
24" (540)	14" - 19.5"
25" (559)	14.5" - 20"
26" (590)	15" - 20"

REAR SEAT HEIGHT d Measure from floor to top of seat tube at back post.

☐ 13" ☐ 13.5" ☐ 14" ☐ 14.5" ☐ 15" ☐ 15.5" ☐ 16" ☐ 16.5" ☐ 17" ☐ 17.5" ☐ 18" ☐ 18.5" ☐ 19" ☐ 19.5" ☐ 20"**STD**

CENTER OF GRAVITY e Measure from front of back post at the seat tube to the center of the rear axle.

☐ 1" ☐ 1.5" ☐ 2" ☐ 2.5" ☐ 3" ☐ 3.5" ☐ 4"**STD**

☐ **Amputee Axle Plate** (Center of Gravity positioned 1" behind back canes. N/A with User Friendly Anti-tips. Limitations will apply)**\$200**

FRONT HANGERS

Composite Footrest with Std Fork					
	3"	4"	5"	6"	8"
85°	•	•			
80°	•	•	•	•	
70°	•	•	•	•	•

Depth Adj Footrest with Std Fork					
	3"	4"	5"	6"	8"
85°	•	•	•	•	
80°	•	•	•	•	
70°	•	•	•	•	•

*Seat to Footrest must be at least 4" less than Front Seat Height

Composite Footrest with Glide Fork					
	3"	4"	5"	6"	8"
85°					
80°					
70°	•	•	•	•	

Depth Adj Footrest with Glide Fork					
	3"	4"	5"	6"	8"
85°					
80°	•	•	•		
70°	•	•	•	•	

☐ 70° ☐ 80° ☐ 85°**STD**

☐ Omit Hangers**NCO**

SEAT TO FOOTREST f Measure from front edge of seat sling to top rear of footrest. Measurement must be at least 2" shorter than Front Seat Height c.

Footrest Type

☐ Composite Flip-up Footrest**STD**

☐ 10.5"- 14" (High Mount Hanger) Set at**\$218**

☐ 14"- 18.5" (Standard Hanger) Set at**NCO**

☐ **Depth and Angle Adjustable Footrest (7"-8.5" will receive 3" extension tube – 9"-10" will receive 5" extension tube)****\$300**

☐ 7"- 10" (Clamp-on High Mount) Set at**\$218**

☐ 10.5"- 14" (High Mount Hanger) Set at**\$218**

☐ 14"- 18.5" (Standard Hanger) Set at**NCO**

☐ Extra-long footplate (9" long)**\$204**

Elevating Leg Rest

Elevating Leg Rest, Single (When choosing a left or right ELR, please select the opposite hanger above. N/A with 14" Seat Width A)

Install on ☐ Left **OR** ☐ Right**\$270**

☐ Elevating Leg Rest, Pair (N/A with 14" Seat Width A)**\$459**

Seat to Footrest/ Elevating Leg Rests f Measure from front edge of seat sling to top rear of ELR footrest.

☐ 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ☐ 19" ☐ 20" ☐ 21"

INTEGRATED RESIDUAL LIMB SUPPORT Mounts in place of front hangers.

- ☐ Residual Limb Support.....\$485
 ☐ 8" x 10" pad **OR** ☐ 8" x 14" pad Install on ☐ L **OR** ☐ R
☐ Residual Limb Support Pair.....\$970
 ☐ 8" x 10" pad **OR** ☐ 8" x 14" pad.

FOOTREST ACCESSORIES

- ☐ Safe Tread Non-Skid Adhesive Patch.....\$22
☐ Heel Loops.....\$60

CALF STRAPS

- Velcro Adjustable.....\$50
 ☐ Large (Recommended for 14"- 17" seat widths.)
 ☐ Extra Large (Recommended for 16"- 22" seat widths.)
 Bodypoint® Padded Velcro Adjustable.....\$117
 ☐ Medium (Recommended for 16"- 20" seat widths.)
 ☐ Large (Recommended for 20"- 22" seat widths.)
 Buckle Adjustable
 ☐ Black Nylon Webbing\$60

SEAT BACK CONFIGURATIONS

SEAT BACK TYPE

- ☐ Fixed Angle, Adjustable Height Seat Back.....STD
☐ Fixed Angle, Adjustable Height Seat Back w/ 8° Bend.....NCO
☐ Folding Adjustable-Angle and Height Seat Back (N/A with HD).....\$295

☐ Folding Stabilizer Bar.....\$457

PUSH HANDLES

- ☐ None.....NCO
☐ Integrated Push Handles.....NCO
☐ Integrated Fold-Down Push Handles (N/A with HD).....\$149
☐ Bolt-on height adjustable Push Handles (N/A with HD).....\$449

SEAT BACK HEIGHT G Measure from top of seat tube at rear of frame to top of backrest upholstery.

Fixed Angle, Adjustable Height Seat Back

Push Handles	Short	Tall	Set at: (Use 1" increments)
None	<input type="checkbox"/> 12.5"- 15.5"	<input type="checkbox"/> 16.5"- 20.5"	Set at: _____
Integrated	<input type="checkbox"/> 12"- 15"	<input type="checkbox"/> 16"- 20"	Set at: _____
Fold-Down	<input type="checkbox"/> 12"- 15"	<input type="checkbox"/> 16"- 20"	Set at: _____
Bolt-on Height Adj.	<input type="checkbox"/> 12.5"- 15.5"	<input type="checkbox"/> 16.5"- 20.5"	Set at: _____

Folding Angle, Adjustable Height Seat Back

Push Handles	Short	Medium	Tall	Set at: (Use .5" increments.)
None	12.5"- 14.5"	15"- 17"	17.5"- 20"	Set at: _____
Integrated	14"- 15.5"	16"- 17"	17.5"- 20"	Set at: _____
Fold-Down	14.5"-15.5"	16"- 17.5"	18"- 20.5"	Set at: _____
Bolt-on Height Adj.	12.5" - 14.5"	15" – 17"	17.5" – 20"	Set at: _____

BACK UPHOLSTERY

- ☐ Standard Back Upholstery (Tension Adjustable. Includes Privacy Flap).....STD
☐ AIR Breathable Back Upholstery (Tension Adjustable. Includes Privacy Flap).....\$385
 ☐ Removable Lumbar Contour Pack ☐ Removable Lateral Contour Pack

☐ Omit Back Upholstery.....NCO

PERMOBIL BACK SUPPORTS (Please see information on page 8 regarding product shipment)

COMFORT ACTA-BACK® Back Support (Backs over 20" wide or tall subject to increase pricing)

- Shell Height ☐ 10" ☐ 12" ☐ 14" ☐ 16" ☐ 18" ☐ 20" ☐ 22".....\$905
 ☐ Comfort-Tek Cover ☐ Stretch-Air Cover.....NCO
 ☐ Comfort Privacy Flap.....\$80

COMFORT ACTA-BACK® Deep Back Support (Backs over 20" wide or tall subject to increase pricing)

Shell Height <input type="checkbox"/> 10" <input type="checkbox"/> 12" <input type="checkbox"/> 14" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/> 22".....	\$905
<input type="checkbox"/> Comfort-Tek Cover <input type="checkbox"/> Stretch-Air Cover.....	NCO
<input type="checkbox"/> 3" Deep Wing <input type="checkbox"/> 5" Deep Wing.....	NCO
<input type="checkbox"/> Comfort Privacy Flap.....	\$80

COMFORT ACTA-EMBRACE® (Backs over 20" wide or tall subject to increase pricing)

Shell Height <input type="checkbox"/> 10" <input type="checkbox"/> 12" <input type="checkbox"/> 14" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/> 22".....	\$1,254
<input type="checkbox"/> Comfort-Tek Cover <input type="checkbox"/> Stretch-Air Cover.....	NCO
<input type="checkbox"/> Comfort Privacy Flap.....	\$80

COMFORT ACTA-RELIEF

Shell Height <input type="checkbox"/> 10" <input type="checkbox"/> 12" <input type="checkbox"/> 14" <input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/> 22".....	\$1,221
<input type="checkbox"/> Comfort-Tek Cover <input type="checkbox"/> Stretch-Air Cover.....	NCO
<input type="checkbox"/> Comfort Privacy Flap.....	\$80

SEAT UPHOLSTERY AND SEATING

<input type="checkbox"/> Tension Adjustable.....	STD
<input type="checkbox"/> Omit Seat Upholstery.....	NCO

PERMOBIL SEAT CUSHIONS (Please see information on page 8 regarding product shipment)

Comfort Cushions (Cushions over 20" wide subject to increase pricing)

<input type="checkbox"/> Comfort Curve Comfort-Tek Cushion.....	\$177
<input type="checkbox"/> Comfort Elements® Comfort-Tek Cushion (Odd size cushions subject to increased pricing).....	\$121
<input type="checkbox"/> Comfort Acta-Embrace™ Comfort-Tek Zero Elevation Cushion.....	\$508
<input type="checkbox"/> Comfort M2 Zero Elevation Cushion with Comfort-Tek cover.....	\$607

ROHO® Cushions (Cushions over 20" wide subject to increase pricing)

<input type="checkbox"/> ROHO Single Compartment <input type="checkbox"/> Low Profile <input type="checkbox"/> Mid Profile <input type="checkbox"/> High Profile	\$713
<input type="checkbox"/> ROHO Quadro Select <input type="checkbox"/> Low Profile <input type="checkbox"/> Mid Profile <input type="checkbox"/> High Profile	\$777
<input type="checkbox"/> ROHO Hybrid Elite <input type="checkbox"/> Single Compartment <input type="checkbox"/> Dual Compartment.....	\$914
<input type="checkbox"/> ROHO Hybrid Select.....	\$1,106

CUSTOM CONFIGURATION OPTIONS

FRONT WHEELS (Please reference table in front hanger section for compatibility)

3" Casters

<input type="checkbox"/> 3"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/ Soft Roll Tire <input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub.....	\$245
--	-------

4" Casters

<input type="checkbox"/> 4"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire.....	STD
<input type="checkbox"/> 4"x.75" Plastic Wheel w/Poly Tire.....	\$70
<input type="checkbox"/> 4"x.75" Light-Up Micro-Caster.....	\$149
<input type="checkbox"/> 4"x1.5" LiteSpeed Billet Aluminum Wheel w/Soft Roll Tire	
<input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub.....	\$259

5" Casters

<input type="checkbox"/> 5"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire.....	STD
<input type="checkbox"/> 5"x1" Plastic Wheel w/Poly Tire.....	\$70
<input type="checkbox"/> 5"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire	
<input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub.....	\$274

6" Casters

<input type="checkbox"/> 6"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire.....	STD
<input type="checkbox"/> 6"x1" Plastic Wheel w/Poly Tire.....	\$70
<input type="checkbox"/> 6"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire	
<input type="checkbox"/> Silver Hub OR <input type="checkbox"/> Black Hub.....	\$274

8" Front Wheels

<input type="checkbox"/> 8"x1" Plastic Wheel w/Poly Tire	NCO
--	-----

FRONT FORKS

<input type="checkbox"/> TiLite Standard Fork	STD
<input type="checkbox"/> Glide Suspension Fork(Max user weight 265 lbs).....	\$596

TiLite X Canadian Order Form
FOLDING FRAME

Effective June 2025 (New items in red)

Page 5 of 8

Rider's Weight: ☐ 125 lbs. ☐ 125-180 lbs. ☐ 180-265 lbs.

REAR WHEELS (See Rear Seat Height Ranges Chart for Rear Seat Height D limitations)

ONE ARM DRIVE

- ☐ One Arm Drive (Requires 24" Rear Wheel, 2" of rear wheel spacing, and 0° Camber, N/A with HD.) (Comes with 24" Spoked rear wheels, MBL SpeedLite (Pneumatic) tires, 0° camber, and Aluminum Silver Ano handrims)..... **\$1,402**
Install on ☐ Left **OR** ☐ Right

STANDARD REAR WHEELS (See Rear Seat Height Ranges Chart for Rear Seat Height D limitations)

	22" (501)	24" (540)	25" (559)	26" (590)
TiLite Wheels	<input type="checkbox"/> STD	<input type="checkbox"/> STD	<input type="checkbox"/> STD	<input type="checkbox"/> STD

When ordering 24" wheel with treaded pneumatic or treaded tire with airless inserts and aluminum black or silver anodized handrims and long tabs, the chair will ship with TiLite standard rear wheel. Otherwise, the chair will ship TiLite Shadow rear wheel.

ALTERNATIVE REAR WHEELS (See Rear Seat Height Ranges Chart for Rear Seat Height D limitations.)

	22" (501)	24" (540)	25" (559)	26" (590)
X-Core Mag®	<input type="checkbox"/> \$269	<input type="checkbox"/> \$269		
Spinergy Spox <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> White	<input type="checkbox"/> \$1,427	<input type="checkbox"/> \$1,427	<input type="checkbox"/> \$1,427	<input type="checkbox"/> \$1,427
Spinergy LX <input type="checkbox"/> Black <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> White <input type="checkbox"/> Bright Orange <input type="checkbox"/> Flamingo Pink <input type="checkbox"/> Purple	<input type="checkbox"/> \$1,547	<input type="checkbox"/> \$1,547	<input type="checkbox"/> \$1,547	<input type="checkbox"/> \$1,547
Spinergy Design Lab LX (Available with black spokes only.) <input type="checkbox"/> Gold Anodize <input type="checkbox"/> Green Anodize <input type="checkbox"/> Blue Anodize <input type="checkbox"/> Red Anodize <input type="checkbox"/> Purple Anodize <input type="checkbox"/> Orange Anodize	<input type="checkbox"/> \$2,382	<input type="checkbox"/> \$2,382	<input type="checkbox"/> \$2,382	<input type="checkbox"/> \$2,382

X-Core Mags requires Aluminum Silver Anodized Handrims with long tabs. N/A with Heavy Duty Frame, Quad Release Axles, Spoke Guards

PNEUMATIC REAR TIRES

- ☐ Treaded pneumatic (N/A with 25"(559) rear wheels)..... **STD**
☐ Treaded with Airless Inserts (N/A with 25"(559) rear wheels)..... **\$157**
☐ MBL SpeedLite..... **STD**
☐ Primo Sentinel..... **\$175**
☐ Schwalbe® Marathon Plus Evolution (N/A with Mag or 22" (501) wheels)..... **\$202**
☐ Primo C-763 Knobby (Available on 24"(540) rear wheels)..... **\$154**

SOLID REAR TIRES

- ☐ Primo Active (Available with 24"(540) and 25"(559) Rear Wheels)..... **\$191**
☐ **Primo Orion (Replaces Primo Express) (N/A with 26"(590) Rear Wheels)**..... **NCO**

CAMBER (Rear wheel spacing will be set to minimum value allowed by your configuration)

- ☐ 0° ☐ 2° ☐ 4°..... **STD**

AXLES

- ☐ Stainless Quick Release..... **STD**
☐ Quad Release (N/A with Mag, 24"(540)/ 25" (559) Spinergy Spox or Spinergy LX Rear Wheels)..... **\$175**
☐ **Fixed Non-removable**..... **\$193**

HANDRIMS

- ☐ Aluminum – ☐ Silver Ano **OR** ☐ Black Ano..... **STD**
☐ Plastic Coated – Black **\$180**
☐ Omit Handrims..... **NCO**

The Natural-Fit® – Original, Larger Profile

- ☐ Standard Grip (Requires Long Tabs)..... **\$454**
☐ Super Grip (Requires Long Tabs)..... **\$514**
☐ Without Thumb Grip..... **\$415**

The Natural-Fit® LT – Reduced Weight, Smaller Profile

- ☐ Standard Grip (Requires Long Tabs)..... **\$546**
☐ Super Grip (Requires Long Tabs)..... **\$585**
☐ Without Thumb Grip..... **\$503**

The Surge (All-In-One Oval with Gription™ Strip)

- ☐ The Surge (Larger Profile)..... **\$591**
☐ The Surge LT (Smaller Profile)..... **\$503**

Q Grip (High-Friction, Anti-Slip Coating)

☐ Q-Grip\$591

TAB LENGTH

☐ Long Tabs.....STD

☐ Short Tabs (N/A with Mag Rear Wheels).....\$47

SPOKE GUARDS

☐ Clear Lexan (N/A with Mag Rear Wheels.).....\$214

WHEEL LOCKS

Aluminum Wheel Locks

☐ Push to LockSTD

☐ Pull to LockSTD

☐ Extension Handles for Aluminum Push/Pull to Lock.....\$80

Composite Wheel Locks

☐ Compact Composite Push to Lock.....\$139

☐ Compact Composite Pull to Lock\$139

☐ Extension Handles for Compact Composite Push/Pull to Lock.....\$111

☐ Compact Composite Scissor Lock.....\$225

Additional Wheel Locks

☐ Grade Aids, Push to Lock (Requires pneumatic tires).....\$472

☐ Hemi Wheel Lock – ☐ Left **OR** ☐ Right.....\$356

SIDE GUARDS

ALUMINUM SIDE GUARDS

Rigid Removable ☐ Small ☐ Standard.....\$318

Quad Removable w/ Hole ☐ Small ☐ Standard.....\$345

ARMRESTS

☐ Omit.....NCO

Standard Armrest

☐ Fixed Height Non-Removable Armrest (N/A with Heavy Duty).....STD

Swing Away Tubular Armrests (N/A HD, Bolt on push handles or Transit/Transport option)

☐ Height Adjustable\$284

Range: ☐ 9" to 12" **OR** ☐ 12" to 15"

☐ Height Adjustable with Sideguards (Side guard will be NCO) (N/A with Heavy Duty).....\$459

☐ 9" to 12" **OR** ☐ 12" to 15"

T-Style Armrests (Requires minimum seat depth of 15")

2-Point Flip Back Height Adjustable T-Style Armrest (N/A with Folding Angle, Adjustable Height Seat Back).....\$410

☐ 8" to 12" **OR** ☐ 10" to 14"

☐ Desk Length Pad (10") **OR** ☐ Full Length Pad (14")

Removeable Height Adjustable T-Style Armrest\$390

☐ 8" to 12" **OR** ☐ 10" to 14"

☐ Desk Length Pad (10") **OR** ☐ Full Length Pad (14")

POSITIONING BELTS

☐ Velcro Adjustable, Black Nylon Webbing.....\$61

☐ Bodypoint® 2" Wide Push-Button Buckle (Auto-Style).....\$66

☐ Bodypoint® 1.5" Wide Push-Button Buckle Padded (Auto-Style, N/A with 19" -20" Rear Widths).....\$147

☐ Bodypoint® 1" Wide Push-Button Buckle Padded (Auto-Style, N/A with 19" -20" Rear Widths).....\$112

☐ Bodypoint® EvoFlex® Pelvic Stabilizer Push Button.....\$174

☐ Q'Straint® WC-19 Wheelchair-Anchored Pelvic Belt (Recommended for use with Transit Tie Down).....\$259

☐ Bodypoint® Optional Band Clamp Mounting Hardware.....\$43

ACCESSORIES

ANTI-TIPS

- ☐ Standard Rear Aluminum - Black.....\$180
☐ User Friendly Flip-up.....\$394

FUNCTIONAL ACCESSORIES

- ☐ Upper Extremity Support Tray
☐ 12" to 14" OR ☐ 13" to 15" Seat Widths A (Requires T Style Armrests.).....\$376
☐ Luggage Carriers (N/A with High Mount Footrests.).....\$189
☐ Crutch Holder - ☐ Left OR ☐ Right.....\$291
☐ Oxygen Tank Holder - ☐ Left OR ☐ Right.....\$356

SOFT GOODS ACCESSORIES

- ☐ Permobil Medical Necessities Bag.....\$123
☐ Permobil Medical Necessities Seat Pouch.....\$79
☐ Bodypoint® Mobility Bag.....\$141
☐ Neoprene Impact Guards, Black.....\$233

GENERAL

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. HCPCS Codes are for reference purposes only. Actual coding and coverage varies by payer. Proper HCPCS coding is not a guarantee of payment. Please verify the appropriate payer's coverage policy.

Get the latest order forms at permobil.com. Prices, standard features, options and accessories are subject to change without notice. All rights reserved. Trademarks are identified by symbols™ and ®. TiLite, TiFit, the TiLite logo, Ti, the Ti logo, and TiShaft are registered trademarks of TiSport, LLC and are licensed to Permobil. Permobil® is a registered trademark of Permobil. Adobe®, Acrobat®, Acrobat Professional® and Reader® are registered trademarks of Adobe Systems Inc. Bodypoint® is a registered trademark of Bodypoint Inc. Kenda® is a registered trademark of Kenda Rubber Industrial Co. Ltd. ROHO®, AGILITY®, Low Profile® and Quadro Select® are registered trademarks of ROHO Inc. Schwalbe® is a registered trademark of Ralf Bohle GmbH. The Natural-Fit® and Out-Front® are registered trademarks of Three Rivers Holdings, LLC. The ACTA-BACK®, ACTA-EMBRACE™, COMFORT-TEK™, COMPASS®, CURVE®, ELEMENTS®, and STRETCH AIR™ are registered trademarks of THE COMFORT COMPANY® and are licensed to Permobil. Q'Straint® is a registered trademark of Q'Straint Inc.


QUOTE THIS CHAIR

If you are using Adobe Acrobat, this order form may be converted into a quote request or into an order by clicking one of the buttons (this feature is available exclusively to Permobil dealers and distributors). A Quote or an Order Confirmation will be sent to you within 24 hours.

ORDER THIS CHAIR

Measurement Guide

- A Seat Width - Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.
B Seat Depth - Measure from front of back post to the front edge of seat sling.
C Front Seat Height - Measure from floor to top of seat tube at front edge of seat sling.
D Rear Seat Height - Measure from floor to top of seat tube at back post.
E Center of Gravity - Measure from front of back post at the seat tube to the center of the rear axle.
F Seat to Footrest - Measure from front edge of seat sling to top rear of footrest
G Seat Back Height - Measure from top of back post to top of seat tube at rear of frame



Seating and Positioning Shipping Information

The following Seating and Positioning products will ship with the chair. Cushions are stocked with Comfort Tek Cover and backs are stocked with Stretch Air Cover. All other Seating and Positioning products will be drop-shipped from the manufacturer. (Width x Depth)

Curve Cushion	Elements Cushion	Acta Embrace Cushion	M2 Cushion Zero Elevation	ROHO Single High Profile Cushion	ROHO Quadtro Cushion
18" x 18"	18" x 18"	18" x 18"	18" x 18"	18" x 18"	18" x 18"
16" x 16"	16" x 18"	16" x 18"	18" x 16"	16" x 18"	16" x 16"
16" x 18"	18" x 16"	16" x 16"	20" x 16"	16" x 16"	20" x 16"
			20" x 20"		20" x 18"
			16" x 16"		

Acta Back	Acta Embrace Back	Acta Relief Back	ROHO Agility Mid Active	ROHO Agility Min Active
18" x 12"	18" x 16"	18" x 16"	18" x 10"	18" x 10"
18" x 16"	16" x 16"	16" x 16"	16" x 10"	16" x 10"
16" x 16"		18" x 12"		
16" x 12"		16" x 12"		