

KNEE CADDY Series



SPECIFICATIONS

25 lbs
17" - 20" Knee pad height
30.5" Width x 21.5" Length
38.5" Handle height
350 lbs weight capacity



FOR THOSE WHO REQUIRE SUPPORT AT THE KNEE, THE KNEE CADDY IS A VERSATILE KNEE SCOOTER WITH A THICK, HEIGHT-ADJUSTABLE KNEE PAD.

Wide-diameter curved steel tubing provides superior stability. Dual rear brake pads actuated by a single brake lever allows for quick and reliable braking..

FEATURES

KNEE PAD Soft, waterproof, durable polyurethane knee pad

COLLAPSIBLE TILLER Folds down easily for quick storage

REMOVABLE BASKET Rigid basket for personal items

HANDLES Ergonomic hand grips

WHEELS 8"x1.5" Polyurethane tires with dual bearings

COLORS Blue



*Net Weight without accessories, backstrap and folding basket



Evolution
TECHNOLOGIES

Toll Free 1 800 556 2558

info@evolution.ca

www.evolution.ca