



- closed plastic arm cuff made of high-quality material
- Rubber tips made of natural rubber with metal inlay
- Special light metal tubes, powder coated in the **COLOR** version in a **COLOUR** of your choice
- 3-fold in the upper part, 13-fold height-adjustable in the lower part with **CLIP**

GRIP OPTIONS

ERGO HARD GRIP

332CEF
Grip length 11 cm



ERGO SOFT GRIP

332CESF
Grip length 11 cm



ERGO WOODEN GRIP

332CEHF
Grip length 12 cm
untreated beechwood



ANATOMIC HARD GRIP

332CEAF
Grip length 12 cm
axially adjustable



ANATOMIC SOFT GRIP

332CEASF
Grip length 12 cm



ANATOMIC WOODEN GRIP

332CEAHF
Grip length 12 cm
untreated beechwood
axially adjustable



The plastic parts are only available in black.
COLOUR SELECTION TUBES

- | | |
|-----------|----------------|
| aubergine | green |
| purple | metallic blue |
| orange | metallic grey |
| pink | metallic red |
| black | textured black |

**FITTING
ACCESSORIES
STARTING ON
PAGE 164**



A 22,5 cm – 27,5 cm
B 72 cm – 102 cm
ø-middle tube 23 mm
ø-lower tube 20 mm

