BRUNO®



ELITE

Curved Outdoor Stair Lift

Made In USA


BRUNO
MADE IN USA



WEATHER-RESISTANT

DURABILITY

Enjoy more time outdoors with the Bruno Elite, the top-selling curved outdoor stair lift in North America.

Handcrafted to precisely match every angle of your stairway, the Bruno Elite curved outdoor stair lift is equipped with marine-grade materials to endure harsh weather and deliver an exceptionally smooth ride.



HIGH-PERFORMANCE RELIABILITY

Weather-Resistant Reliability
Engineered to withstand extreme climates and temperatures (0˚F - 125˚F, -18˚C - 52˚C), the Elite is designed with marine-grade materials to deliver year-round performance.

Handcrafted Quality
The Elite's custom-made rail delivers a precise fit and smooth ride, while the covered gear rack helps shield the drive system from the elements.

Dependable
Battery-powered design delivers top-to-bottom operation even in the event of a power outage.

Saves Space
Vertical rail maximizes open space on steps. Armrests, seat and footrest fold flat when not in use.

Easy to Use
Armrest control glides the stair lift up and down. Includes two wireless remotes.





AVAILABLE OPTIONS

Add extended rail positions to park the stair lift away from the steps or out of view, or an optional mid-park charging station for stairs with middle landings.

Opt for a wider seat or larger foot-rest to further customize your Elite stair lift chair.

Not all options shown.







★★★★★

"I love my Bruno stair lift. It works so smoothly. **It is truly a blessing.**"

Sharon W.
Bruno Stair Lift Customer







TECH SPECS

Lift Capacity:
400 lb (181 kg)

Standard Seat:
18 in wide x 15 in deep (457 mm x 381 mm)

Larger Seat Option:
22.5 in wide x 16.5 in deep
(572 mm) x (429 mm)

Standard Footrest:
10.5 in long x 9.4 in wide
(260 mm x 232 mm)

Larger Footrest Option:
12 in long x 11 in wide (305 mm x 279 mm)

Vertical Rail:
Installs as close as 7.5 in (191 mm) to the wall or edge of staircase

WARRANTY

Limited 5-Year Warranty on Major Components (motor, gearbox, drivetrain and stationary rail system) for the lifetime of the Bruno stair lift as long as the original purchaser owns the stair lift.

2-Year Warranty on all parts after date of installation by an Authorized Bruno Dealer.

30-Days Labor Coverage after date of installation by an Authorized Bruno Dealer for any manufacturing defect.

ABOUT BRUNO

Veteran-Founded | **Family-Owned** | **Made in USA**

Bruno's mission is to help make life easier for people with mobility challenges by delivering the highest quality stair lifts, scooter lifts, platform lifts and home elevators.

BRUNO®

© Bruno Independent Living Aids, Inc. 2025
2025-09-09-elite-outdoor-curve-mkt-sl-5.indd